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2000006355

KAN TORBASI 3LÜ T&amp;T CPDA-1 450 ML

**TANIM**

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
## DESCRIPTION

Closed blood bag systems made of PVC, consisting of a sampling system, 16 g needle, transfer bags, platelet storage bag (TOTM) and including CPDA-1 anticoagulant solution.

## INTENDED USE

For separation, collection, storage, transporting and processing of red blood cells (RBCs), plasma and platelets (PLTs) from Platelet Rich Plasma (PRP) per "Guide to the preparation, use and quality assurance of Blood Components".

## INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

## INDICATIONS FOR USE

This product is intended to be used for collection of whole blood and preparing blood components from whole blood for all population to aid any treatment. Blood components can be used for restoring blood components, improving oxygen delivery to tissues, restoring clotting factors and reducing or preventing bleeding according to the treatment decision of physician.

## CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

## PRECAUTIONS

- Do not use if needle protector is opened.
- Do not add medication to blood.
- Cross match the contents of the bag with recipient before transfusion.
- Use transfusion sets in compliance with EN ISO 1135-4.
- Follow validated blood component preparation procedures per local requirements.
- The tubes are suitable to be connected with a Sterile Connecting Device.
- It is intended for human blood and components.
- During shipment, up to two shipping boxes may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them, and, they should be loaded in such a way to prevent them from falling over.
- Single use only, after usage it should be thrown as a Medical Waste. **Caution:** Risk of infection when used more than once.
- For the explanation of symbols, see Symbols Description Table.

**Caution:** This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

## GENERAL PREPARATIONS BEFORE USAGE

- Open the outer wrap at the tear notches and remove one blood bag system, check for any visible defects of tubes, needle, anticoagulant and additive solutions in terms of appropriate volume, discoloration or particulate contaminants.

**Caution:** Do not use defective sets, make sure the solutions are clear.

- Apply blood pressure cuff or tourniquet on donor's arm, identify venipuncture site and release the cuff or tourniquet.
- Disinfect the venipuncture site. Cover the area with sterilized gauze strip and do not touch until the time of venipuncture.

**Caution:** Beware of iodine sensitive donors!

- Signs on the needle protector and hub should be the in the same line (Figure 9). Apply some pressure to determine any leakages.
- Place the collection bag below the level of donor's arm.

**Caution:** While opening the blood bag system, blood taking tube should be held from needle protector not from the needle safety cover. Pay attention not to lock the needle safety cover before collection of blood by mistake.

## BLOOD COLLECTION (indicative)

- Before removing the needle's cap, close the clamp on the donor line connected to collection bag (for systems in Figure 1 and 2). If the bag system has a donor line breakaway part, there is no need to close the clamp (Figure 3 and 4).
- Reapply tourniquet or inflate blood pressure cuff until previously selected vein is prominent again.
- Hold the needle from marked side. To remove the needle cap, twist to break the seal between the needle and its cap and remove it by sliding it down the axis of the needle. **Caution:** Be sure not to touch the cannula!
- Apply the venipuncture immediately. Venipuncture is very important for a non-coagulated blood.
- Strap the needle to the donor's arm.


Figure 9


Figure 1

Figure 2

Figure 3

Figure 4

Figure 5

6. Follow the following steps according to the sampling bag system:

**a. System with Sampling Welder**

- i Break the breakaway part on the tube to the sampling holder (for system Figure 1, No: 2).
- ii Insert the sampling tubes into the tube holder and collect the required blood sample (Figure 1, No: 3).  
**Caution:** Use bigger tubes first.
- iii After the sampling is completed, open the clamp on the tubing to the bag and start blood collection (Figure 1, No: 1).
- iv Follow step 7.

**b. System with Sampling Bag**

- i Take initial blood to the sampling bag (Figure 2, No: 5).
- ii Close the clamp on the tubing to the sampling bag (Figure 2-3-4, No: 2).
- iii Break the breakaway part on the tubing to the bag by turning left and right (Figure 7) and blood taking starts.
- iv Place the sampling bag into the palm tube holder facing down. Insert the sampling tubes into the tube holder connected to the sampling bag and collect the required blood sample. **Caution:** Use bigger tubes first.  
**Caution:** Blood sample in the sampling bag should immediately transferred to the tubes, otherwise they may get coagulated.
- v After sampling complete, break the breakaway part on the tubing to the collection bag (Figure 3 and 4, No: 1). Follow step 7.

Figure 6

7. Hang the tube with the tube holder onto donor's wrist in a way not to compress the needle.

8. Mix blood and anticoagulant gently and periodically (approximately in every 45 seconds) during collection. If a mixer or a collection system is used, please refer to their operator's manual.

9. Make sure there is continuous blood flow. If there is continuous, adequate blood flow and constant agitation, blood collection is expected to be completed within 12 minutes. **Caution:** If it is above 12 minutes, collected blood may not be suitable for Platelets, Fresh Frozen Plasma or Cryoprecipitate AHF.

10. Collect the quantity of blood within the limits indicated on the bag label. Monitor the blood being collected.

11. After blood collection is completed, stop the blood mixing device and loosen the tourniquet or deflate the pressure cuff.

12. For the blood bag systems with needle safety cover, while holding the needle with one hand hold the needle safety cover with the other one and then leave the needle, pull the needle to the safety cover by pulling the tube under the safety cover and be sure that it is locked (Figure 5, No: 6). **Caution:** To prevent accidental injuries or contamination of the personnel, insert the needle into the locked system or throw it to the medical waste.

13. Close the blood taking tube under the Y part by tube sealing device and remove the needle from the system. **Caution:** Throw locked needle kit to the medical waste.

14. Strip donor tubing as much as possible through the bag and start sealing. Work quickly to prevent the blood from clotting in the tubing. Invert bag several times to mix thoroughly; then allow tubing to refill with anticoagulated blood from the bag. Repeat this procedure twice (Figure 6).

15. Check the blood bag again for any faults.

16. Make sure all donor information required is filled out on the label on the blood bag.

17. After collection, whole blood should be placed at controlled room temperature (T) ( $20 < T < 24$  °C) for a min of 2 hours and a max of 6 hours with 45° angle, if more or overnight holding period is used it can be stored at 1-6 °C in a blood bank refrigerator. For gaining fresh frozen plasma and platelets, complete all the steps within 8 hours after collection. In order to reduce red blood cell contamination, hold the bag upright by the ports and flick the ports thoroughly with the fingers until all traces of RBC are cleared back into the bag.


Figure 7

**COMPONENT SEPARATION (indicative)**

1. Load the whole blood filled bag system in the centrifuge buckets. **Caution:** Proper positioning of the bags will prevent burst of the bags during centrifugation.

2. Centrifuge at 2500 g for 15 minutes under 22 °C. **Caution:** Centrifuge settings may vary according to the centrifuge brand used.

3. Place the centrifuged bag into the blood component extractor.

4. After the centrifugation, place the system into the plasma separator.

5. Break the break-off part at the opening of the bag by turning left and right (Figure 7).

6. Squeeze the tubes with clamps except TOTM bag to be used for plasma.

7. Transfer plasma to the TOTM bag marked with "5" via an extractor. Perform the plasma separation as described in the extractor instruction manual.

8. Close the tube of red blood cell bag with welder.

9. Store the red blood cells immediately as specified at storage section.

10. Centrifuge plasma at 3000 g for 10 minutes again. **Caution:** Centrifuge settings may vary according to the centrifuge brand used.

11. After centrifugation, PLTs remain at the bottom of the bag (shaped like an arch) and platelet poor plasma on the top.

12. Place TOTM bag in the extractor and unclamp the plasma bag.

13. Transfer the platelet poor plasma on the top to the empty plasma bag.

14. Close the tube of platelet bag with welder and store as specified in the storage part.

15. Close the tube of platelet poor plasma bag with welder and store as specified in the storage part.

16. Freeze the platelet poor plasma at -40 °C.

Figure 8

**STORAGE**

1. RBCs should be placed into blood bank refrigerators designed for this process under 1-6 °C.

2. Store the plasma in the freezer. If the separation takes place within 6 hours after collection, it can be stored at -18 °C/-25 °C for 3 months, -30 °C/-40 °C for 6 months and at lower temperatures for 12 months.

3. If PLTs are transferred into the special PVC (TOTM) bag, they can be stored up to 5 days under incubator agitation at 22 °C.

**TRANSFUSION**

1. Inspect the container and blood for any defects.

2. Before blood transfusion, follow the validated local procedure.

3. Before transfusion, make sure that serological test results are acceptable and cross match test results with the recipient. If appropriate, apply cross match again at bedside.

4. Platelets should be inspected before transfusion to ensure that no platelet aggregates are visible.

5. Hold the top part of the T part with one hand and bag from the bottom part of T part with other hand, and twist of the upper part by 270° (Figure 8).

6. Keep holding the bag from T part, remove the cap of the transfusion spike with other hand.

7. Insert the spike by ¾ of its length into outlet port. Do not touch the spike of the transfusion set and open outlet port to prevent contamination.

8. After inserting the spike, squeeze the bag to the level that transfusion set filter is filled with blood.

9. Hang the bag onto the stand. The unit is ready for transfusion.

SYMBOL	TÜRKÇE	ENGLISH	FRANCAIS	ITALIANO	DEUTSCH	POLSKI	ESPAÑOL	PORTUGUÊS	DANSK	SVENSKA	ΕΛΛΗΝΙΚΑ	BULGARIAN	SLOVENSKY
	XX mL Kan ya da Kan Bileşeni Kabı + nominal hacim	Blood or blood component container + nominal volume	Recipient contenant du sang ou des composants sanguins + volume nominal	Sacca sangue o emponente + volume nominale	Behälter für Blut oder Blutkomponenten + Nennvolumen	Pojemnik na krew lub składniki krwi + objętość nominalna	Bolsa de sangre o componente + volumen nominal	Saco da sangue ou componente + volume nominal	Blod- eller blodkomponentbeholder + nominelt volum	Blodkomponent-påse + nominelt volum	Περιέκτης αίματος ή συστατικού αίματος + ονομαστικός όγκος	Контейнер за човешка крв или крвини съставки + номинален обем	Vak na krv alebo zložku krvi + nominálny objem
	Trombosit	Platelets	Plaquettes	Piastrine	Thrombozyten	Płytki krwi	Plaquetas	Plaquetas	Trombocytter	Trombocytter	Λιπιδεμένα	Тромбоцити	Krvné doštičky
	Tam kan bileşeni kabı	Whole blood	Sang total	Sangue intero	Vollblut	Krew	Sangre total	Sangue total	Fuld blod	Helblod	Περιέκτης ολικού αίματος	Цяла кръв	Celá krv
	Kırmızı kan hücreleri konsantresi	Red blood cell concentrate	Erythrocytes	Emazie concentrată	Erythrozyten	Erytrocyty	Hematies	Eritrócitos	Erytrocytter	Erytrocyter	Συμπυκνωμα Ερυθροκυττάρων	Еритроцитен концентрат	Koncentrát červených krviniek
	Buffy Coat	Buffy coat	Couche leucoplaquettaire	Buffy Coat	Buffy-Coat Beutel	Kozuszek leukocytarno-płytkowym	Capa leucoplaquetária	Camada leucoplaquetária	Buffyccoat	Buffy Coat	Buffy coat	Крватен серум	Buffy-coat
	Trombosit kabı	Platelets container	Poche de conservation du concentré plaquettaire	Sacca per piastrine	Thrombozyten Beutel	Pojemnik z krvinekami płytkowymi	Bolsa para Plaquetas	Recipiente de plaquetas	Trombocytbeværningspose	Trombocyt påse	Περιέκτης αιμοπετάλιων	Сак за съхранение на тромбоцитен концентрат	Vak na krvné doštičky
	Kırmızı kan hücreleri kabı	Red blood cell container	Poche contenant les Erythrocytes	Sacca per la conservazione delle emazie concentrate	Erythrozyten Beutel	Pojemnik czerwonych krwinek	Bolsa para Hematies	Recipiente de ematies	Pose til Erytrocytter	Erytrocyter påse	Περιέκτης ερυθροκυττάρων	Сак за съхранение на еритроцитен концентрат	Vak na červené krvinky
	Plazma kabı	Plasma container	Poche contenant le plasma	Sacca per la conservazione del plasma	Plasma Beutel	Plasma container	Bolsa para Plasma	Recipiente de plasma	Pose til plasma	Plasma påse	Περιέκτης πλάσματος	Сак за съхранение на плазма	Vak na plazmu
	Buffy Coat kabı	Buffy coat container	Poche contenant la couche leucoplaquettaire	Sacca per buffy coat	Buffy-Coat Beutel	Pojemnik z kozuszkami leukocytarno-płytkowym	Bolsa para "Buffy coat"	Recipiente de Camada leucoplaquetária	Pose til buffyccoat	Buffy coat påse	Περιέκτης Buffy coat	Сак за съхранение на крватен серум	Vak na buffy-coat
	Tam kan bileşeni kabı	Whole blood container	Poche de sang total	Sacca per sangue intero	Vollblutbeutel	Caly worek krwi	Contenedor de sangre total	Recipiente de sangue total	Pose til fuld blod	Helblodspåse	Περιέκτης ολικού αίματος	Сак за съхранение на цяла кръв	Vak na celú krv
	İşlem kabı	Processing container	Poche de transfert	Sacca di trasferimento	Transfer Beutel	Pojemnik do przetwarzania	Bolsa de transferencia	Recipiente de transferencia	Transførspose	Överföringspåse	Περιέκτης επεξεργασίας	Контейнер за преработка	Vak na spracovanie
	Antikoagulan solüsyon	Anticoagulant solution	Solution anticoagulante	Soluzione anticoagulante	Antikoagulans-Lösung	Rozwiązanie przeciwzakrzepowe	Solución anticoagulante	Solução anticoagulante	Antikoagulerende løsning	Antikoagulanlösning	αντιπηκτικό διάλυμα	Антикоагуляционен разтвор	Antikoagulačný roztok
	Katkü çözeltisi	Additive solution	Solution additive	Soluzione conservante	Additivlösung	Rozwiązanie dodatki	Solución aditiva	Solução aditiva	Additiv løsning	Tillsättslösning	Προσθετο διάλυμα	Добавен разтвор	Prídavný roztok
	Lökosit filtreme	Leucocyte filtration	Filtration des leucocytes	Filtrazione leucocitaria	Leukozytenfiltration	Filtracja leukocytów	Filtración de leucocitos	Filtração de leucócitos	Leukocytfiltrering	Leukocytfiltering	Διήθηση των Λευκοκυττάρων	Филтрация на левкоцити	Filtrácia leukocytov
	Diş ambalajda göze görülebilen bozuma belirtilir varsa kullanmayınız	Do not use if there is visible sign of deterioration	Ne pas utiliser si des signes de détérioration sont visibles	Non utilizzare in caso di segni visibili di deterioramento	Nicht verwenden, wenn Anzeichen des Verfalls sichtbar sind	Nie używać, jeśli nosisi ślady uszkodzenia	No utilizar si hay signos visibles de deterioro	Não utilizar se forem visíveis sinais de deterioração	Må ikke anvendes, hvis der er synlige tegn på nedbrydning	Använd ej vid synliga tecken på försämring	Μην το χρησιμοποιείτε εάν υπάρχει ορατή ένδειξη αλλοίωσης	Да не се използва, ако има видими белези на увреждане	Nepoužívať pri viditeľných známkach poškodenja
	Doğal kauçuk lateks içermez	Free of natural rubber latex	Dépourvu de latex naturel	Privo di lattice naturale	Latexfrei	Nie zawiera naturalnego lateksu	Libre de látex de caucho natural	Não contém látex de borracha natural	Frit for naturgummilætex	Latexfri	Δεν περιέχει φυσικό λάτεξ	Не съдържа латекс	Neobsahuje prírodný kaučukový latex
	DEHP içerir ya da bannılır.	Contains or presence of DEHP	Trace ou présence de phthalates : phthalate de di-2-éthylhexyle (DEHP)	Contiene ftalati o presenta tracce di ftalati Di-(2-etilil)ftalato (DEHP)	Enthält DEHP	Zawiera ftalan	Contiene o puede que haya ftalatos presentes Di-(2-ethyl)ftalato (DEHP)	Presença ou conteúdo de DEHP	Produkt eller emballage indeholder DEHP	Innehåller eller förekommer av DEHP	Περιέχει φυσικό DEHP ή ftalyl	Съдържа или наличие на DEHP	Obsahuje alebo môže byť prítomný DEHP
	Hava vermeyiniz.	Do not vent	Ne pas utiliser de prise d'air.	Non introdurre aria	Nicht belüften	Przechowywać szczelnie zamknięte	No ventilar	Não ventilar	Må ikke udluftes	Får ej ventileres	Μην χρησιμοποιείτε αεραγωγό	Да не се пробива	Neodvzdovávať
	Radyasyon ile svi yolu sterilizasyon	Sterile fluid path by radiation	Trajet stérile. Stérilisé par irradiation.	Percurso del liquido sterile. Sterilizzato per irradiazione.	Sterile Fluidweg durch Strahlung	Nie stosować, jeśli opakowanie jest uszkodzone	Via de fluido estéril por RADIACION	Percurso do fluido esterilizado por radiação	Steril vaskbænde ved bestråling	Steril vâtskebana genom strålning	Αποστειρωμένη διαδρομή υγρού από ακτινοβολία	Стерилен поток чрез радиация	Dráha tekutiny je sterilizovaná žiarením
	Isi ile svi yolu sterilizasyon	Sterile fluid path by STEAM	Trajet de fluide stérilisé par VAPEUR	Percurso del liquido sterile. Sterilizzato per Ossido di etilene	Sterile Fluidweg durch Dampf	Sterylnie droga przepływu pary przez parę	Via de fluido estéril por Vapor	Percurso do fluido esterilizado pelo vapor	Steril vaskbænde med damp	Steril vâtskebana genom ÅNGA	Αποστειρωμένη διαδρομή υγρού από τον ατμό	Стерилен поток чрез парна стерилизация	Dráha tekutiny je sterilizovaná parou
	Hasar görmüş yitirmiş veya delinmiş paketeri kullanmayınız.	Do not use if package is damaged	Ne pas utiliser si le produit, sa barrière de stérilité ou son emballage a été endommagé ou s'il présente des signes d'altération.	Non utilizzare se il prodotto, la sua barriera sterile o il suo contenimentamento sono danneggiati o mostrano segni di deterioramento	Nicht verwenden, wenn die Verpackung beschädigt ist	Sterylnie droga przepływu pary przez parę	No usar si el producto, su barrera estéril o su envase está dañado o muestra cualquier signo de deterioro	Må ikke anvendes, hvis produktet, sterilbarrieren eller pakningerne er beskadigede eller der er tegn på brud	Använd inte om förpackningen är skadad	Μην το χρησιμοποιείτε εάν η συσκευασία έχει υποστεί φθορά	Не използвайте, ако опаковката е с нарушавана целост	Nepoužívať, ak je obal poškodený	
	Bu sıcaklık arası saklayınız.	Store between these temperature range	Stockez dans cette plage de température	Conservare in questo intervallo di tempo	Aufbewahren in diesem Temperaturbereich	Przechowywać w tym zakresie temperatur	Almacenar dentro de este rango de temperatura	Armazenar dentro deste intervalo de temperatura	Opbevares ved dette temperaturområde	Lagras inom detta temperaturområde	Αποθηκεύεται μεταξύ της συγκεκριμένης θερμοκρασίας	Да се съхранява при следните температури	Skladujte v tomto teplotnom rozsahu
	Alüminyum paket açılması takiben X gün içinde kullanılmadır.	Contents of Al pouch must be used within X days of opening	Le contenu de l'emballage en Aluminium doit être utilisé dans un délai de X jours après ouverture	Il contenuto della confezione di Al va utilizzato entro X giorni dall'apertura	Der Inhalt der Verpackung muss innerhalb von X Tagen nach dem Öffnen verwendet werden	Zawartość Al opakowania musi być użyta w ciągu X dni od momentu otwarcia opakowania	El contenido de envases de aluminio debe utilizarse en los X días siguientes a su apertura	O conteúdo da embalagem de Al tem de ser utilizado no prazo de X dias após a sua abertura	Når Al-pakningen er åbnet, skal indholdet anvendes inden for X dage	Innehållet i Al-förpackningen måste användas inom X dagar efter öppnandet	Το περιεχόμενο της συσκευασίας αλουμινίου πρέπει να χρησιμοποιηθεί εντός Χ ημερών από το άνοιγμα	Съдържанието на опаковката да се използва до Х дни след отваряне	Obsah hliníkového fólieového vaku sa musí spotrebovať do X dní od otvorenia.
	Paketi içeriği açıldıktan sonra 10 gün içinde kullanılmalıdır.	Contents of pack must be used within 10 days of opening	Le contenu de l'emballage doit être utilisé dans un délai de 10 jours après ouverture	Il contenuto della confezione di Al va utilizzato entro 10 giorni dall'apertura	Der Inhalt der Verpackung muss innerhalb von 10 Tagen nach dem Öffnen verwendet werden	Zawartość opakowania musi być użyta w ciągu 10 dni od momentu otwarcia opakowania	El contenido del envase debe utilizarse en los 10 días siguientes a su apertura	O conteúdo do envase deve ser utilizado no prazo de 10 dias após a sua abertura	Når pakningen er åbnet, skal indholdet anvendes inden for 10 dage	Innehållet i förpackningen måste användas inom 10 dagar efter öppnandet	Το περιεχόμενο της συσκευασίας πρέπει να χρησιμοποιηθεί εντός 10 ημερών από το άνοιγμα	Съдържанието на опаковката да се използва до 10 дни след отваряне	Obsah vrečka sa musí spotrebovať do 10 dní od otvorenia
	Aprojen svi yolu	Non-pyrogenic fluid path	Trajet apyrogène	Percurso del liquido apirogeno	Nichtpyrogenen Flüssigkeitspfad	Apirogena ścieżka przepływu	Paso de fluido aprotégeno	Percurso de fluido aprotégeno	Non pyrogen vaskbænde	Icke-pyrogen vâtskebana	Μη πυρογόνο διαδρομή υγρού	Непирогенен поток	Dráha tekutiny je nepyrogenou
	Alet veya paket	pcs per package.	unités par emballage	pz per confezione	Stück pro Packung	Szt. w opakowaniu	Unidades que contiene el paquete	peças por pacote	stk pr pakke	st per förpackning	τεμάχια ανά ή συσκευασία	Брой в опаковка	ks v balení.
	Depolama işleminde 8 kolitten fazla üst üste konulmamalıdır.	Stacking limit by number. Maximum number of identical packages which may be stacked on one another, where "8" is the limiting number	Ne pas empiler plus de 8 cartons	Sovraporre al max 8 cartoni/scatole	Stapeln Sie nicht mehr als 8 Boxen	Dopuszczalna ilość warstw piętrzenia: "8"	Límite máximo de estiva 8 cajas	Não empilhar mais de 8 caixas	Må ikke stables mere end 8 kasser	Stapla inte mer än 8 boxar	Μην στοιβάζετε περισσότερο από 8 κιβώτια κατά την αποθήκευση	Складвайте в зависимост от броя продукти. Максимален брой пакети, които могат да бъдат подредени един върху друг, когато "8" е определят максимален брой	Krajný počet pri stávaní. Maximálny počet rovnakých balení, ktoré sa môžu stávať jedno na druhé, pričom "8" je krajná hodnota
	Tek kullanılmıktır.	Do not re-use.	Ne pas réutiliser	Monouso	Nicht wiederverwenden	Nie używać ponownie.	No reutilizar	Não reutilizar	Må ikke genbruges	Ateranvänd inte.	Πρόν μιας χρήσης	Да не се използва повторно	Nevhodné na opätovné použitie.
	Kullanmadan önce kullanma kılavuzunu okuyunuz.	Read instructions for use.	Consulter la notice	Consultare le istruzioni per l'uso	Lesen Sie die Gebrauchsanweisung für die Anwendung	Przeczytaj instrukcję użycia	Leer las instrucciones de uso	Leia as instruções de utilização	Læs brugsanvisningen	Läs bruksanvisningen	Διαβάστε τις οδηγίες χρήσης πριν τη χρήση	Запознайте се с инструкцията за употреба	Prečítajte si pokyny na používanie.
	Dikkat.	Caution, read accompanying documents.	Attention, consulter les documents d'accompagnement.	Attenzione, consultare i documenti all'interno	Vorsicht, lesen Sie die Begleitdokumentation	Uwaga, przeczytaj dołączoną dokumentację	Atención, leer los documentos adjuntos	Cuidado, leia os documentos de acompanhamento	Se brugsanvisningen	Varning, läs medföljande dokument	Προσοχή, διαβάστε τα συνοδευτικά έγγραφα	Внимание, прочетете придружаващите документи	Pozor, prečítajte si sprievodné dokumenty.
	Direkt güneş ışığı ve yüksek sıcaklığı maruz bırakmayınız.	Keep away from sunlight. Keep away from heat.	Protéger des rayons directs du soleil et d'une source puissante de chaleur.	Proteggere dalla luce solare diretta. Mantenere lontano da fonti di calore	Vor Sonnenlicht schützen. Hitze fernhalten	Chronić przed nagrzaniem (cieplem)	Proteger del sol. Mantener lejos del calor	Mantênhã longe da luz solar. Manter afastado do calor	Beskyt mod direkte sollys og varme	Hålla borta från solljus. Hålla borta från värme	Να φυλάσσεται μακριά από το φως του ήλιου. Να φυλάσσεται μακριά από θερμότητα	Да не се излага на слънчева светлина. Да не се сержва	Chránit pred priamým svetlom. Chránit pred teplotou.
	Kuru yerde saklayınız. Name maruz bırakmayınız.	Keep dry.	Conservar à l'abri de la pluie	Conservare in luogo asciutto	Trocken lagern	Chronić przed wilgocią	Conservar en lugar seco	Manter seco	Opbevares tørt	Förvaras torrt	Διατηρήσει το στεγνό	Да се съхранява на сухо място	Skladujte v suchu.
	Kolaylıkla kırılabilir, dikkatli tutunuz.	Fragile, handle with care	Fragile, manipuler avec précautions.	Facilmente divisibile, tenere ben ferma.	Verfallsdatum	Snadno rozbídné, držte ji opatrně	Fragil. Manipule con cuidado	Facilmente quebrável, segure-o com cuidado	Let i stykker, hold det forsigtigt	Utgångsdatum	Ευθραστο, χειριστείτε με προσοχή	Чупливо, внимавайте при употреба	Krehké, zaobchádzajte opatrne
	Son kullanma tarihi.	Use by, date.	Date d'expiration	Scadenza	Vor scharfen Gegenständen entfernt halten	Nie używać po	Fecha de vencimiento	Data de validade	Utløbsdato	Hålla borta från vassa föremål	Ημερομηνία λήξης	Използвайте преди	Spotřebujte do dátumu.
	Torbaları kesici ve delici materyallerden uzak tutunuz	Keep away from sharp objects.	Tenir à l'écart des objets tranchants	Tenere lontano da oggetti appuntiti o taglienti	Vor scharfen Gegenständen entfernt halten	Trychmy z dala od ostrzych przedmiotów	Mantener alejado de objetos afilados	Manter longe de objetos cortantes	Holdes væk fra skarpe genstande	Förvaringspås för trombocytter upp till fem dagar	Να φυλάσσεται μακριά από αιχμηρά αντικείμενα	Παζете от остри предмети	Chránit pred ostrými predmetmi.
	5 güne kadar Trombosit saklama torbasi	PLT storage bag up to five days	Poche de conservation de plaquettes jusqu'à cinq jours	Sacca per la conservazione delle piastrine fino a cinque giorni	Aufbewahrungsbeutel für Thrombozyten bis zu fünf Tage	Pojemnik do przechowywania płytek do pięciu dni	Bolsa para el almacenamiento de plaquetas hasta cinco días	Recipiente de armazenamento de plaquetas até cinco dias	Opbevaringspose for trombocytter op til fem dage	Förvaringspås för trombocytter upp till fem dagar	Περιέκτης αποθήκευσης αιμοπετάλιων έως πέντε ημέρες	Сак за съхранение на тромбоцитен концентрат до пет дни	Vak na uschovu doštickek do petih dni
	7 güne kadar Trombosit saklama torbasi	PLT storage bag up to seven days	Poche de conservation de plaquettes jusqu'à sept jours	Sacca per la conservazione delle piastrine fino a sette giorni	Aufbewahrungsbeutel für Thrombozyten bis zu sieben Tage	Pojemnik do przechowywania płytek do siedmiu dni	Bolsa para el almacenamiento de plaquetas hasta siete días	Recipiente de armazenamento de plaquetas até sete dias	Opbevaringspose for trombocytter op til syv dage	Förvaringspås för trombocytter upp till sju dagar	Περιέκτης αποθήκευσης αιμοπετάλιων έως επτά ημέρες	Сак за съхранение на тромбоцитен концентрат до седем дни	Vak na uschovu doštickek do sedmich dni
	Numune alma bölgisi	Sampling site	Site d'échantillonnage	Siti di campionamento	Probenahemstelle	Punkt poboru próbek	Para toma de muestras	Local de amostragem	Prøvetagning sted	Provtagningsställe	Σημείο δειγματοληψίας	сайт за вземане на проби	Miesto odberu
	Tekrar sterilize etmeyiniz.	Do not re-sterilize	Ne pas restériliser	Non ristilizzare	Nicht erneut sterilisieren	Nie należy ponownie sterylizować	No reesterilizar	Não reesterilize	Må ikke reesteriliseres	Resteriliseras inte	Μην αποστειρώνετε ξανά	Не реставрирайте	Nesmie sa opätovne sterilizovať
	Üretici	Manufacturer	Fabricant	fabricante	Hersteller	Producent	Fabricante	Fabricante	Fabrikant	Tillverkare	Κατασκευαστής	Производител	Výrobca
	Seri Numarası	Batch Code	Code du lot	Codice batch	Batch-Code	Kod partii	Código de lote	Código de lote	Batchcode	Batchkod	Κωδικός Παρτίδας	Код на партидата	Číslo šarže
	Ürün Kodu	Catalogue Number	Numéro de catalogue	Numero di catalogo	Katalognummer	Numer katalogowy	Número de catalogo	Catálogo de número	Katalognummer	Katalognummer	Αριθμός Καταλόγου	Каталожен номер	Katologové číslo
	CE Uygulama İlaresi ve Belgelendirme Kuruluş Numarası	CE Marking of Conformity and Identification Number of Notified Body	Marquage CE de conformité et numéro d'identification de l'organisme notifié	Marchio CE di conformità e numero d'identificazione dell'organismo notificato	CE-Konformitätsbescheinigung und Kennnummer der benannten Stelle	Oznakowanie CE zgodności i numeru identyfikacyjnego jednostki notyfikowanej	Marcação CE de conformidade e Número de Identificação do Organismo Notificado	Marcação CE de conformidade e Número de identificação do Organismo Notificado	CE-mærkning af overensstemmelse og identifikation Antal benyttede organer	CE-märkning av överensstämmelse och identifikation Antal benyttade organ	Επίσημο CE συμμόρφωσης και αριθμός αναγνωριστικού του κοινοποιημένου οργανισμού	СЕ маркировка за съответствие и идентификационен номер на нотифициран орган	Označenie zhody CE a identifikačné číslo notifikačného orgánu