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2*****6732	KAN TORBASI T&B ERİTROSİT LÖKOSİT FİLTRELİ DÖRTLÜ CPD/SAGM 450 ML NAS2
22204406745	KAN TORBASI T&B ERİTROSİT LÖKOSİT FİLTRELİ BEŞLİ CPD/SAGM 450 ML NAS2
22004006745	KAN TORBASI T&B ERİTROSİT LÖKOSİT FİLTRELİ BEŞLİ CPD/SAGM 450 ML
22104406745	KAN TORBASI T&B ERİTROSİT LÖKOSİT FİLTRELİ CPD/SAGM 450 ML

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

Şekil 9

### APPLICABLE CODES


2*****6732	BLOOD BAG T&B WITH RED BLOOD CELL LEUKOCYTE FILTER QUADRUPLE CPD/SAGM 450 ML PDS2
22204406745	BLOOD BAG T&B WITH RED BLOOD CELL LEUKOCYTE FILTER QUINTUPLE CPD/SAGM 450 ML PDS2
22004006745	BLOOD BAG T&B WITH RED BLOOD CELL LEUKOCYTE FILTER QUINTUPLE CPD/SAGM 450 ML
22104406745	BLOOD BAG T&B WITH RED BLOOD CELL LEUKOCYTE FILTER CPD/SAGM 450 ML

\* means applicable for all systems with same last four digit


### DESCRIPTION

A blood bag system made of PVC, including CPD anticoagulant solution and SAG-M additive solution, connected in a closed system with integrated red cell (RC) leukocyte filter, blood sampling system, 16 g needle, transfer bags, platelet (PLT) storage bag (TOTM).

#### A. 2\*\*\*\*\*6732

#### B. 2\*\*\*\*\*6745

### INTENDED USE

For collection, separation and storage of leukodepleted red cells, plasma and platelets per "Guide to the preparation, use and quality assurance of Blood Components."

### PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

### INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

### INDICATIONS FOR USE

This product is intended to be used for collection of whole blood and preparing leukocyte-reduced blood components from whole blood for all population to aid any treatment. Blood components can be used for restoring blood components, improving oxygen delivery to tissues, restoring clotting factors and reducing or preventing bleeding according to the treatment decision of physician.

### CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

### PRECAUTIONS

1. Do not use if needle protector is opened.
2. Do not add medication to blood.
3. Cross match the contents of the bag with recipient before transfusion.
4. Use transfusion sets in compliance with EN ISO 1135-4.
5. Follow validated blood component preparation procedures per local requirements.
6. The tubes are suitable to be connected with a Sterile Connecting Device.
7. It is intended for human blood and components.
8. During shipment, up to two shipping boxes may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them, and they should be loaded in such a way to prevent them from falling over.
9. Single use only, after usage it should be thrown as a Medical Waste. **Caution:** Risk of infection when used more than once.
10. For the explanation of symbols, see Symbols Description Table.

**Caution:** This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

### GENERAL PREPARATIONS BEFORE USAGE

1. Open the outer wrap at the tear notches and remove one blood bag system, check for any visible defects of tubes, needle, anticoagulant and additive solutions in terms of appropriate volume, discoloration or particulate contaminants.

**Caution:** Do not use defective sets, make sure the solutions are clear.

**Caution:** After sticking the barcode label with the donor information on the bag, apply pressure by hand to make sure the label is adhered.

2. Apply blood pressure cuff or tourniquet on donor's arm, identify venipuncture site and release the cuff or tourniquet.
3. Disinfect the venipuncture site. Cover the area with sterilized gauze strip and do not touch until the time of venipuncture.

**Caution:** Beware of iodine sensitive donors!

4. Signs on the needle protector and hub should be in the same line (Figure 11). Apply some pressure to determine any leakages.
5. Place the collection bag below the level of donor's arm.

**Caution:** While opening the blood bag system, blood taking tube should be held from needle protector not from the needle safety cover. Pay attention not to lock the needle safety cover before collection of blood by mistake.

### BLOOD COLLECTION (indicative)

1. Before removing the needle's cap, close the clamp on the donor line connected to collection bag (for systems in Figure 1 and 2). If the bag system has a donor line breakaway part there is no need to close the clamp (Figure 3 and 4).
2. Reapply tourniquet or inflate blood pressure cuff until previously selected vein is prominent again.
3. Hold the needle from marked side. To remove the needle cap, twist to break the seal between the needle and its cap and remove it by sliding it down the axis of the needle. **Caution:** Be sure not to touch the cannula!
4. Apply the venipuncture immediately. Venipuncture is very important for a non-coagulated blood.
5. Strap the needle to the donor's arm.
6. Follow the following steps according to the sampling bag system:


Figure 11


Figure 1

Figure 2

Figure 3

Figure 4

#### a. System with Sampling Tube

- i Break the breakaway part on the tube to the sampling holder (for system Figure 1, No: 2).
- ii Insert the sampling tubes into the tube holder and collect the required blood sample (Figure 1, No: 3).  
**Caution:** Use bigger tubes first.
- iii After sampling complete, open the clamp on the tubing to the bag and blood collection starts (Figure 1, No: 1).
- iv Follow step 7.

#### b. System with Sampling Bag

- i Take initial blood to the sampling bag (Figure 2, No: 5),
- ii Close the clamp on the tubing to the sampling bag (Figure 2-3-4, No: 2),
- iii Break the breakaway part on the tubing (Figure 3, 4) to the bag by turning left and right (Figure 7) and blood taking starts.
- iv Place the sampling bag in the palm tube holder facing down. Insert the sampling tubes into tube holder connected to the sampling bag and collect the required blood sample. **Caution:** Use bigger tubes first.  
**Caution:** Sample blood in the sampling bag should immediately transferred to the tubes otherwise they may get coagulated.
- v After sampling complete, break the breakaway part on the tubing to the collection bag (Figure 3 and 4, No:1).
- vi Follow step 7.


Figure 5

7. Hang the tube with the tube holder onto donor's wrist in a way not to compress the needle.
8. Mix blood and anticoagulant gently and periodically (approximately in every 45 seconds) during collection. If a mixer or a collection system is used, please refer to their operator's manual.
9. Make sure there is continuous blood flow. If there is continuous, adequate blood flow and constant agitation, blood collection is expected to be completed under 12 minutes.  
**Caution:** If it is above 12 minutes, collected blood may not be suitable for Platelets, Fresh Frozen Plasma or Cryoprecipitate AHF.
10. Collect the quantity of blood within the limits indicated on the bag label. Monitor the blood being collected.
11. After blood collection is completed, stop the blood mixing device and loosen the tourniquet or deflate the pressure cuff.
12. For the blood bag systems with needle safety cover, while holding the needle with one hand hold the needle safety cover with the other one, then leave the needle, pull the needle to the safety cover by pulling the tube under the safety cover and be sure that it is locked (Figure 5, No:6). **Caution:** To prevent accidental injuries or contamination of the personnel, insert the needle into the locked system or throw it to the medical waste.
13. Close the blood taking tube from under the Y part with tube sealing device and remove the needle from the system. **Caution:** Throw locked needle kit to the medical waste.
14. Strip donor tubing as much as possible through bag, start sealing. Work quickly to prevent the blood from clotting in the tubing. Invert bag several times to mix thoroughly; then allow tubing to refill with anticoagulated blood from the bag. Repeat this procedure twice (Figure 6).
15. Check the blood bag again for any faults.
16. Make sure all the donor information required is filled out on the label of the blood bag.
17. After collection, whole blood should be placed at controlled room temperature (T) ( $20^{\circ}\text{C} < \text{T} < 24^{\circ}\text{C}$ ) for a min of 2 hours and a max of 6 hours with  $45^{\circ}$  angle, if more or overnight holding period is used it can be stored at  $-1-6^{\circ}\text{C}$  in a blood bank refrigerator. For gaining fresh frozen plasma and platelets, finish all the steps within 8 hours after collection.

### COMPONENT SEPARATION (indicative)

#### Top&Bottom Quadruple&Quintuple Inline Bag Systems (Intended for Plasma and PLT separation by Buffy Coat (BC))

1. Place the bags into centrifuge bucket with the breakaway part stands on above.  
**Caution:** Proper positioning of the bags will prevent burst of the bags during centrifugation.
2. If the blood will be used as a whole blood, break the breakaway part on the main bag and transfer whole blood into 500 ml capacity transfer bag.
3. Centrifuge at 3500 g for 15 min at  $22^{\circ}\text{C}$ . **Caution:** Centrifuge settings may vary according to the centrifuge brand used.
4. Place the centrifuged bag to the blood component extractor. If quintuple bag is used, clamp the PLT bag's tube.
5. Use the blood component separator according to its separator using manual.
6. Firstly, break the main bag's breakaway part, then the transfer bags (Figure 7).
7. Separator transfers plasma into the plasma bag, red cells to the transfer bag while keeping BC in the main bag.
8. When the procedure is done, seal transfer bags of plasma and red cells, and remove them from the system.
9. Store plasma as identified at the point 14.
10. If quintuple bag is used, lock main bag with BC with PLT clamp, rinse the bag and rest it for an hour.
11. Centrifuge main bag and platelet (5) transfer bag at 500 g for 6 min at  $20^{\circ}\text{C}$ . **Caution:** Centrifuge settings may vary according to the centrifuge brand used.
12. Separator transfer PLTs to the transfer bag. **Caution:** PLTs are separated after centrifuge.
13. Platelets can be stored at  $22^{\circ}\text{C}$  under agitation in an incubator up to 5 days in PLT bag made of TOTM. **Caution:** Place the bags horizontally and do not put on top of each other.
14. Store the plasma in the deep freeze. If the separation takes place within 6 hours after collection, it can be stored at  $-18^{\circ}\text{C}/-25^{\circ}\text{C}$  for 3 months,  $-30^{\circ}\text{C}/-40^{\circ}\text{C}$  for 6 months and at lower temperatures for 12 months.

Figure 6

### FILTRATION

1. Hook the SAG-M bag to filter hanger from 150 cm high (Figure 8).
2. Break the breakaway part of SAG-M bag by turning it left and right (Figure 7).
3. SAG-M solution flows through the filter into the RCs bag.
4. After SAG-M process is finished lock the tube under leukocyte filter with clamp.  
**Caution:** Beware to put the clamp to the closest point of RCs bag.
5. Mix the RCs bag gently before hook up. **Caution:** Proper mixing the bag decreases filtration time.
6. Reverse the whole kit, hook RCs to the top and leave the SAG-M bag on the bottom (Figure 9).
7. Open the clamp. RCs flow into the filter and leukocytes are filtered and they are collected in the SAG-M bag.
8. Filtration time is 25 min average at room temperature.
9. Invert bag several times to mix thoroughly; then allow tubing to refill with RCs from the bag. Repeat this procedure twice (Figure 6).

### TRANSFUSION

1. Inspect the container and blood for any defects.
2. Before blood transfusion, follow the validated local procedures.
3. Before transfusion, make sure that serological test results are acceptable, cross match test results match with the recipient. If appropriate, apply cross match again at bedside.
4. Platelets should be inspected before transfusion to ensure that no platelet aggregates are visible.
5. Hold the top part of the T part with one hand and bag from the bottom part of T part with other hand, and twist of the upper part by  $270^{\circ}$  (Figure 10).
6. Keep holding the bag from T part, remove the cap of the transfusion spike with other hand.
7. Insert the spike by  $\frac{3}{4}$  of its length into outlet port. Do not touch the spike of the transfusion set and open outlet port to prevent contamination.
8. After inserting the spike, squeeze the bag to the level that transfusion set filter is filled with blood.
9. Hang the bag from its holder. The unit is ready for transfusion


Figure 7


Figure 8

Figure 9

Figure 10

KANSUK LABORATUARI SANAYI VE TICARET A.Ş.  
Yassıören Mahallesi Fırat Sk. No:14/1 Arnavutköy-İstanbul/Turquie  
Tél: +90 212 592 15 76, Fax: +90 212 580 37 72

[www.kansuk.com](http://www.kansuk.com), E-mail: [kansuk@kansuk.com](mailto:kansuk@kansuk.com)

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