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2*****6384	KAN TORBASI DÖRTLÜ CPD/SAGM 450 ML
2*****6361	KAN TORBASI ÜÇLÜ CPD/SAGM 450 ML

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
*means applicable for all systems with same last four digit


DESCRIPTION

A Blood Bag system made of PVC, including CPD anticoagulant solution and SAG-M additive solution, blood sampling system, 16 G needle, transfer bags, platelet (PLT) storage bag (TOTM).

A. 2*****6384

B. 2*****6361

INTENDED USE

For collection, separation and storage of, Plasma and Platelets (PLT) from Platelet Rich Plasma (PRP) per "Guide to the preparation, use and quality assurance of Blood Components."

PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

INDICATIONS FOR USE

This product is intended to be used for collection of whole blood and preparing blood components from whole blood for all population to aid any treatment. Blood components can be used for restoring blood components, improving oxygen delivery to tissues, restoring clotting factors and reducing or preventing bleeding according to the treatment decision of physician.

CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

PRECAUTIONS

- Do not use if needle protector is opened.
- Do not add medication to blood.
- Cross match the contents of the bag with recipient before transfusion.
- Use transfusion sets in compliance with EN ISO 1135-4.
- Follow validated blood component preparation procedures per local requirements.
- The tubes are suitable to be connected with a Sterile Connecting Device.
- It is intended for human blood and components.
- During shipment, up to two shipping boxes may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them, and they should be loaded in such a way to prevent them from falling over.
- Single use only, after usage it should be thrown as a Medical Waste. **Caution:** Risk of infection when used more than once.
- For the explanation of symbols, see Symbols Description Table.

Caution: This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

GENERAL PREPARATIONS BEFORE USAGE

- Open the outer wrap at the tear notches and remove one blood bag system, check for any visible defects of tubes, needle, anticoagulant and additive solutions in terms of appropriate volume, discoloration or particulate contaminants.

Caution: Do not use defective sets, make sure the solutions are clear.

- Apply blood pressure cuff or tourniquet on donor's arm, identify venipuncture site and release the cuff or tourniquet.
- Disinfect the venipuncture site. Cover the area with sterilized gauze strip and do not touch until the time of venipuncture.

Caution: Beware of iodine sensitive donors!

- Signs on the needle protector and hub should be in the same line (Figure 9). Apply some pressure to determine any leakages.
- Place the collection bag below the level of donor's arm.

Caution: While opening the blood bag system, blood taking tube should be held from needle protector not from the needle safety cover. Pay attention not to lock the needle safety cover before collection of blood by mistake.

BLOOD COLLECTION (indicative)

- Before removing the needle's cap close the clamp on the donor line connected to collection bag (for systems in Figure 1 and 2). If the bag system has the breakaway valve into the donor line, there is no need to close the clamp (Figure 3 and 4).
- Reapply tourniquet or inflate blood pressure cuff until previously selected vein is prominent again.
- Hold the needle from marked side. To remove the needle cap, twist to break the seal between the needle and its cap and remove it by sliding it down the


Figure 9


Figure 1 Figure 2 Figure 3 Figure 4

Figure 5

- axis of the needle. **Caution:** Be sure not to touch the cannula!
4. Apply the venipuncture immediately. Venipuncture is very important for a non-coagulated blood.
 5. Remove the needle to the donor's arm.
 6. Follow the following below steps according to the sampling bag system:

a. System with Sampling Welder

- i Open the breakaway valve on the tube to the sampling holder (for system Figure 1, No: 2).
 - ii Insert the sampling tubes into the tube holder and collect the required blood sample (Figure 1, No: 3).
- Caution:** Use bigger tubes first.
- iii After once the sampling is completed, open the clamp on the tubing to the bag and blood collection starts (Figure 1, No: 1).
 - iv Follow step 7.

Figure 6

b. System with Sampling Bag

- i Take initial blood to the sampling bag (Figure 2, No: 5).
 - ii Close the clamp on the tubing to the sampling bag (Figure 2-3-4, No: 2).
 - iii Break the breakaway part on the tubing to the bag by turning left and right (Figure 7) and blood taking starts.
 - iv Place the sampling bag into the palm tube holder facing down. Insert the sampling tubes into tube holder connected to the sampling bag and collect the required blood sample. **Caution:** Use bigger tubes first.
- Caution:** Blood sample in the sampling bag should immediately transferred to the tubes otherwise they may get coagulated.
- v After sampling complete, open the breakaway valve on the tubing to the collection bag (Figure 3 and 4, No: 1). Follow step 7.


Figure 7

7. Hang the tube with the tube holder onto donor's wrist in a way not to compress the needle.
8. Mix blood and anticoagulant gently and periodically (approximately in every 45 seconds) during collection. If a mixer or a collection system is used, please refer to their operator's manual.
9. Make sure there is continuous blood flow. If there is continuous, adequate blood flow and constant agitation, blood collection is expected to be completed within 12 minutes. **Caution:** If it is above 12 minutes, collected blood may not be suitable for Platelets, Fresh Frozen Plasma or Cryoprecipitate AHF.
10. Collect the quantity of blood within the limits indicated on the bag label. Monitor the blood being collected.
11. After blood collection is completed, stop the blood mixing device and loosen the tourniquet or deflate the pressure cuff.
12. For the blood bag systems with needle safety cover, while holding the needle with one hand hold the needle safety cover with the other one, then leave the needle, pull the needle to the safety cover by pulling the tube under the safety cover and be sure that it is locked (Figure 5, No: 6). **Caution:** To prevent accidental injuries or contamination of the personnel, insert the needle into the locked system or throw it to the medical waste.
13. Close the blood taking tube downstream the Y part by tube sealing device and remove the needle from the system. **Caution:** Throw locked needle kit to the medical waste.
14. Strip donor tubing as much as possible through bag, starting by seal. Work quickly, to prevent the blood from clotting in the tubing. Invert bag several times to mix thoroughly; then allow tubing to refill with anticoagulated blood from the bag. Repeat this procedure twice (Figure 6).
15. Check the blood bag again for any faults.
16. Make sure all the donor information required is filled out on the label of the blood bag.
17. After collection, whole blood should be placed at controlled room temperature (T) ($20 < T < 24$ °C) for a min of 2 hours and a max of 6 hours with 45° angle, if more or overnight holding period is used it can be stored at 1-6 °C in a blood bank refrigerator. For gaining fresh frozen plasma and platelets, finish all the steps within 8 hours after collection. In order to reduce red cell contamination, hold the bag upright by the ports and flick the ports thoroughly with the fingers until all traces of red cells (RCs) are cleared back into the bag.

COMPONENT SEPARATION

CENTRIFUGATION

1. Load the full blood filled bag system in the centrifuge buckets. **Caution:** Proper positioning of the bags will prevent burst of the bags during centrifugation.
2. Centrifuge at 2500 g for 15 minutes under +22 °C 2500 g. **Caution:** Centrifuge settings may vary according to the centrifuge brand used.
3. Place the centrifuged bag into the blood component extractor.

Triple Bag (Used for separation of red cells, plasma and PLT's by PRP)

4. After the centrifugation, carefully place the centrifuged unit into the manual or automated plasma separator.
5. Break the break-off part at the opening of the bag by turning left and right. (Figure 7)
6. Squeeze the tubes with clamps except plasma separator.
7. Transfer the plasma to the transfer bag marked with "5" with the help of extractor. Perform plasma separation as specified in the extractor guide.
8. If platelet separation will not be performed, close the plasma tube with welder and remove from the system. Store the plasma as specified in the storage part.
9. If PLT's will be separated, read component separation part, item number 13.
10. Break the breakaway part of SAG-M bag by turning it left and right. SAG-M solution flows through the filter into the RCs bag.
11. Seal the tube of RBCs bag.
12. Store the red cells immediately as specified in the storage part.
13. Plasma is re-centrifuged at 3000 g for 10 min. **Caution:** Centrifuge settings may vary according to the centrifuge brand used.
14. After centrifuge, PLTs remain at the bottom, and plasma on the top of the bag.
15. Place the plasma bag to the blood component extractor, and unlock the clamp.
16. Transfer plasma to the empty transfer bag.
17. Close the tube of platelets bag with welder and store the platelets as specified in the storage part.
18. Close the tube of the bag of platelet-poor plasma with welder and store as specified in the storage part.

Figure 8

Quadruple Bag (Used for separation of red cells, plasma and PLT's and cryoprecipitate by PRP)


19. Use the same instruction of number 1 to 16 in the component separation part.
20. Store the platelets as specified in the storage part.
21. Freeze plasma purged of platelets under -40 °C.
22. Frozen plasma is then thawed.
23. Centrifuge at 3000 g for 10 minutes. **Caution:** Centrifuge settings may vary according to the centrifuge brand used.
24. After centrifuge supernatant part of the crystalized plasma is transferred into other bag. A jelly part will remain at the bottom which is called cryoprecipitate.


STORAGE

1. RCs should be placed into blood bank refrigerators under 1-6 °C.
2. Store the plasma in the freezer. If the separation takes place within 6 hours after collection, it can be stored at -18 °C / -25 °C for 3 months, -30 °C / -40 °C for 6 months and at lower temperatures for 12 months.
3. If it is transferred into the special PVC (TOTM) bag, PLTs can be stored up to 5 days under incubator agitation at 22 °C.
4. Freeze the cryoprecipitate. The storage date will be the same of the plasma.

TRANSFUSION

1. Inspect the container and blood for any defects.
2. Before transfusing the RCs, make sure that serological test results are acceptable, cross match test results match with the recipient. If appropriate, apply cross match again at bedside.
3. Platelets should be inspected before transfusion to ensure that no platelet aggregates are visible.
4. Hold the top part of the T part with one hand and bag from the bottom part of T part with other hand, and twist of the upper part by 270° (Figure 8).
5. Keep holding the bag from T part, remove the cap of the transfusion spike with other hand.
6. Insert the spike by ¾ of its length into outlet port. Do not touch the spike of the transfusion set and open outlet port to prevent contamination.
7. After inserting the spike, squeeze the bag to the level that transfusion set filter is filled with the relevant blood component.
8. Hang the bag onto the stand. The unit is ready for transfusion.

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