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

DESCRIPTION

Closed blood bag system made of PVC consisting of a sampling system, 16 g needle, bags and including CPD as anticoagulant solution and SAG-M as additive solution.

INTENDED USE

Used for collection, separation and storage of whole blood.

PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

INDICATIONS FOR USE

This product is intended to be used for collection of whole blood and preparing blood components (plasma and platelet) from whole blood for pediatric population to aid any treatment. Blood components can be used for restoring blood components, improving oxygen delivery to tissues, restoring clotting factors according to the treatment decision of physician.

CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

GENERAL PREPARATIONS BEFORE USAGE

1. Open the outer wrap of the blood bag through the tearing notch and remove a blood bag system.
2. Attach the sphygmomanometer or tourniquet on donor's arm, identify venipuncture site and release the tourniquet or sphygmomanometer.
3. Disinfect the venipuncture site. Cover the area with sterilized gauze strip and do not touch until the time of venipuncture.
 - CAUTION:** Beware of iodine sensitive donors!
4. Before entering the vein, review the blood bag system, tubings and needle for visible errors and color change.
5. Signs on the needle protector and hub should be in the same line. Apply some pressure to determine any leakages.
6. Place the collection bag below the level of donor's arm.
7. Check the anticoagulant solution and additional solutions for the correct volume, color and particle content.

CAUTION: While opening the blood bag system, blood taking tube should be held from needle protector not from the needle safety cover. Pay attention not to lock the needle safety cover before collection of blood by mistake.

INSTRUCTION FOR USE

BLOOD COLLECTION:

IF BLOOD BAG SYSTEMS WITH SAMPLING ARE USED, SEE POINTS 1-4; IF THE SYSTEM DOES NOT HAVE SAMPLING, CONTINUE FROM POINT 5.

1. Before removing the needle's cap, close the clamp on the tubing reaching the collection bag. If the system's blood collection line has breakaway valves, it is not necessary to close the clamp off.
 2. Break open the breakaway valve on the tubing to the sampling holder / bag (Figure 1, No: 2; Figure 2, No: 3). **NOTE:** There are not any breakaway valves on the tubing to the sampling tubing holder / bag in the systems with a blood collection line with breakaway valves.
 3. Depending on the sampling system used, follow one of the following processes:
 - c. **If a system with tube holder is used;**
 - i. Insert the sampling tubes into the tube holder and collect the required blood sample (Figure 1, No: 3). 3).
 - ii. **CAUTION:** First, samples are taken into bigger tubes.
 - iii. Open the clamp on the tubing to the bag and start blood collection (Figure 1, No: 1).
 - d. **If a system with sample bags is used;**
 - i. Take the blood to the sampling bag (Figure 2, No: 5).
 - ii. Close the clamp on the tubing to the sampling bag (Figure 2, No: 2).
 - iii. Open the clamp on the tubing to the bag and start blood collection (Figure 2, No: 1).
 - iv. Place the sampling bag in the palm tube holder facing down. Insert the sampling tubes into the tube holder connected to the sampling tube and collect the required blood sample (Figure 2, No: 4). **CAUTION:** Blood sample in the sampling bag should immediately transferred to the tubes; otherwise, the blood may get coagulated.
- NOTE:** After the sampling process is finished, break the breakaway valve on the tubing to the sampling tubing bag in the systems with a blood collection line with breakaway valves.
4. Hang the tubing with the tubing holder onto donor's wrist in a way not to compress the needle. Continue from point 9.
 5. If the system used does not have sampling, attach a temporary clamp or knot the blood collection tubing.
 6. Until the previously selected vein is visible again, inflate the sphygmomanometer again or apply the tourniquet.
 7. Hold the needle from marked side. To remove the needle casing, turn the cap in a circular motion so that the adhesion zone between the needle and the casing is broken and scrape down along the needle axis. **CAUTION:** Be sure not to touch the needle!
 8. Apply the venipuncture with the needle immediately. It is very important to apply the venipuncture cleanly and smooth for a complete blood collection without clots.
 9. Tape the needle to the donor's arm as it will hold a needle.
 10. Mix the blood and anticoagulant solution slowly and periodically (approximately every 45 seconds). If a blood collection device or any collection system is used, follow the instruction manuals.
 11. Make sure there is continuous blood flow. If there is continuous, adequate blood flow and constant agitation, blood collection is expected to be completed within 12 minutes. **CAUTION:** If it is above 12 minutes, collected blood may not be suitable for Platelets, Fresh Frozen Plasma or Cryoprecipitate AHF separation.
 12. Collect the quantity of blood within the limits indicated on the bag label. Monitor the blood being collected.
 13. When the blood collection process is complete, the blood collection shaker takes the blood at the set weight or milliliter; then stop the blood collection by attaching a clamp to the blood collection tubing.
 14. Stop the blood collection shaker device, untie the turnstile after;
 15. The clamp of the blood collection shaker is turned on or lower the pressure of the sphygmomanometer.
 16. When the blood collection process is completed, to remove the needle, hold it with the thumb and index fingers of the other hand while holding the sterile gauze in one hand, and remove the needle completely without applying pressure. Let the donor gently press the sterile gauze to the venipuncture point.
 17. Holding the needle with one hand, wait for the blood from the blood collection tubing to flow into the bag so that the needle tip is upwards.
 18. After this process, close the clamp on the tubing to the bag. In systems without sampling, knot the tubing. **CAUTION:** If the clamp does not close, there may be a risk of blood gushing out of the needle protection cover while the needle is locked.


Figure 1

Figure 2


Figure 3

19. For the blood bag systems with needle safety cover, while holding the needle with one hand hold the needle safety cover with the other one and then leave the needle, pull the needle to the safety cover by pulling the tubing under the safety cover and be sure that it is locked (Figure 3, No: 6). **CAUTION:** To prevent accidental injuries or contamination of the personnel in case of systems without needle safety cover, insert the needle into the needle destruction machine or throw it to the bio-danger box.
20. Close the blood taking tubing under the Y part by tubing sealing device and remove the needle from the system. **CAUTION:** Throw locked needle system to the medical waste.
21. Scrape the blood from the blood collection tubing completely from the place of origin to the bag. Work quickly to prevent the blood from clotting in the tubing. Stir the bag by inverting it over a few times, then allow tubing to refill with anticoagulated blood. Repeat this procedure twice (Figure 4).
22. For cross-comparison purposes, segment the tubing connected to the collection bag with the tubing shut-off machine so that the segment numbers are clearly and fully readable. The segment closest to the collection bag is always attached to the bag for identification purposes.
23. Check the blood bag again for any faults.
24. Make sure that all donor information required is filled out on the label on the blood bag.
25. After complete blood collection, immediately place it in a blood bank cabinet designed for this process at 1-6 °C. If platelets are to be separated, do not cool the blood, but leave it at an angle of 45° until the temperature reaches 20-24 °C (min one hour). For gaining fresh frozen plasma and platelets, complete all the steps within 8 hours after collection.

Figure 4

Quad Pediatric Component Bag (used for Erythrocyte and Plasma separation and Separation of Erythrocytes into low volume bags)

1. Load the whole blood filled bag system in the centrifuge buckets. **Caution:** Proper positioning of the bags will prevent burst of the bags during centrifugation.
2. Centrifuge at 2500 g for 15 minutes under +4 °C. **CAUTION:** Centrifuge settings may vary according to the centrifuge brand used.
3. Place the centrifuged bag into the blood component extractor.
4. Break the break-off part at the opening of the bag by turning left and right (Figure 5).
5. Perform the plasma separation as described in the extractor instruction manual.
6. Transfer the plasma to the transfer bag. Disconnect the tubings from the system by closing them with the welding machine.
7. If the SAG-M bag is used, break the fracture part of the SAG-M bag and add the SAG-M solution to the main bag on top of the erythrocytes.
8. Close the tubing of the erythrocyte bag with welding machine and immediately place the erythrocytes to be kept in a blood bank cabinet designed for this process at 1-6 °C.
9. Insert a temporary knot or clamp on the tubing that connects the erythrocyte bag to the other transfer bags.
10. The entire system should be immediately placed in a blood bank cabinet designed for this purpose to be stored at 1-6°C.
11. Repeat the points 7-13 as the demand for low-volume erythrocytes arrives until the empty bags are depleted.
12. Store the plasma in the freezer. If the separation takes place within 6 hours after collection, it can be stored at -18 °C/-25°C for 3 months, -30°C/-40°C for 6 months and at lower temperatures for 12 months.


Figure 5

BLOOD DONATION

1. Inspect the container and blood for any faults.
2. Before applying the blood, check whether the serological tests applied are appropriate and whether the results of the cross-comparison test are consistent with the blood given. If appropriate, apply cross match again at bedside.
3. Use a transfusion kit in accordance with EN ISO 1135-4.
4. Hold the top part of the T part with one hand and bag from the bottom part of T part in the bag with other hand, and twist of the upper part by 270° (Figure 6).
5. Remove the cap of the transfusion tip with your other hand while holding the lower part of the T part in the bag.
6. Insert the spike by ¾ of its length into outlet port. Do not touch the spike of the transfusion set and open outlet port to prevent contamination.
7. After inserting the spike, squeeze the bag to the level that transfusion set filter is filled with blood. Hang the bag onto the stand from the bottom.

Figure 6

GENERAL WARNINGS

During storage and transportation:

1. During transportation, without throwing the blood bag boxes, place the boxes in the car so that no more than two boxes are stacked together, and, if they are carried by hand or by transport trolley, so that no heavy load will fall on them.
2. In the storage process, stack up to eight boxes so that they do not overlap and do not tip over.

Before opening the outer packaging:


1. Store blood bags with outer packaging in a clean, dry and cool place.
2. Do not use if there is a visible deterioration in the outer packaging.
3. Blood bags are sterile packed.

After the outer packaging is opened:

1. Store the bags in a cool, dry environment inside the outer packaging.
2. Use blood bags within 10 days after the outer packaging is opened.

Before use:

1. Blood bags should be used in accordance with the instructions for use.
2. Intended for human blood and components.
3. Sterile and apyrogen - Sterilized with steam.
4. Do not use bags whose solution is not clear.
5. Do not use if there is a visible deterioration in the bags.
6. Keep the bag and lines away from cutting objects.
7. Do not use if needle protector is opened. The lines on the needle should be directly opposite each other (Figure 7).
8. Do not close after opening the needle cover.
9. Do not use if there are no liquid route covers or if they are loose.
10. Do not aerate.
11. Store the blood between 1-6 °C.
12. Do not add medicine to the blood bag.
13. Compare the blood in the bag and the recipient's blood with the Cross test before the transfusion.
14. Use a transfusion kit in accordance with EN ISO 1135-4.
15. Mix before use and, if necessary, bring to the appropriate temperature under controlled conditions.
16. When frozen, plastic is more fragile.

Figure

After use

After usage it should be thrown as a Medical Waste. **CAUTION:** Risk of infection when used repeatedly.

CAUTION: This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

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