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20013068150

LABORATUVAR YIKANMIŞ ERİTROSİT ( RBC) LÖKOREDÜKSİYON FİLTRE

**TANIM**

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**APPLICABLE PRODUCT**

20013068150	LABORATORY LEUKOREDUCTION FILTER for WASHED RBC
-------------	---

**DESCRIPTION**

They are disposable, sterile transfusion sets which contain a transfer spike, di-(2-ethylhexyl)-phthalate (DEHP) plasticized poly(vinyl chloride) (PVC) tubing, two blood bags, a leukocyte filter, drip chamber and check valve (Figure 1). They are sterilized by Gamma irradiation.


Figure 1

**INTENDED USE**

For preparation and washing of leukocyte-depleted red blood cell (RBC) suspensions with isotonic saline solution and storage according to the "Guide to the Preparation, Use and Quality Assurance of Blood Components".

**PATIENT POPULATION**

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

**INTENDED USER AND INTENDED USE ENVIRONMENT**

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

**INDICATIONS FOR USE**

This product is intended to be used for leukocyte reduction and washing of red blood cells for all population to aid any treatment. Red blood cells can be used for restoring blood components and improving oxygen delivery to tissues according to the treatment decision of physician.

**CONTRAINDICATIONS**

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

## PRECAUTIONS

- The tubes are suitable for connection with a Sterile Connecting Device.
- It is intended to be used for human blood and blood components. During shipment, up to two containers may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them, and, they should be loaded in such a way to prevent them from falling over.
- Single use only. After use, should be disposed of as a Medical Waste. **Caution:** Risk of infection when used more than once.
- For explanation of symbols, see the Symbols Description Table.

**Caution:** This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. However, medical procedures should not be avoided as the benefits outweigh any potential health risks associated with exposure to DEHP. Please refer to the current literature to make an informed decision.

## PROCEDURE STEPS

### A. GENERAL PREPARATIONS BEFORE USE

1. Remove the filter system from the primary packaging. The filter system is ready to be connected to the blood bag (Figure 1).
2. Please visually inspect for any damage to the filter system or if a component is lost or missing. Do not use defective sets.
3. Close the clamps in the tubing after the spike and under the filter (Figure 1, no. 5).
4. Filter should be connected to the bag on a proper bench.

### B. USER INSTRUCTIONS

1. Connect the filter system to the blood bag to be filtered.
  - a. If the connection will be performed with a sterile connecting device, the tube of the RBC unit and tube of the filter (under the spike) (Figure 1, no. 3d) have to be connected. **Caution:** See the user manual of the sterile connecting device.
  - b. If the connection will be performed by spiking, open the safety cap of the filter spike (Figure 1, no. 6), open the safety cap of the outlet port of the blood bag, then insert the spike firmly keeping the bag flat on the bench. (Figure 2) **Caution:** See the user manual of the blood bag.)

Figure 2

2. Hang the blood bag (Figure 3).

Figure 3

3. Open the clamps on the line of the filter (Figure 1, no. 5).
4. Prime the filter (Figure 4).


Figure 4

5. Blood will flow through the filter into the transfer bag.
6. After 90% of the blood is filtered, close the clamp (Figure 1, no. 5).

Figure 5

7. Squeeze the transfer bag and push the air out through the check valve until the RBCs reach the bag, then re-open the clamp and let the RBCs flow through the tube and the filter until the back of the housing is cleared. (Figure 5).
8. When all blood has passed through the filter, remove the filtered RBCs from the system by sealing the tube with a sealer (between no. 10 and no. 15, Figure 1).
9. Insert the spike (Figure 1, no. 8) into the saline container and fill the RBC bag with saline solution.
10. Remove the RBC bag filled with saline solution by sealing the tube (Figure 1, no. 3d) at the closest point to the bag. Centrifuge according to the validated internal procedures.
11. After the centrifugation, place the RBC bag on a plasma extractor, open the clamp on the tubing of the waste bag (Figure 1, no. 5) and transfer the supernatant to the waste bag (Figure 6).


Figure 6

12. Remove the waste bag from the system by sealing the tubing with a sealer (Figure 1, no. 3f) at the closest point to the RBC bag.

### C. STORAGE

Storage time of the washed RBC should be as short as possible and should not exceed 24 hours.

### D. TRANSFUSION

1. Use transfusion sets in compliance with EN ISO 1135-4.
2. Inspect the container and/or blood component for any defects.
3. Before transfusing the blood, make sure that the serological test results are acceptable and the cross-match test results show compatibility with the recipient. If appropriate, repeat the cross-match test at bedside.
4. Hold the upper part of the T part of the bag with one hand and the bottom part with the other, and twist the upper part by 270°.
5. Keep holding the bag from the T part, and remove the cap of the transfusion spike with the other hand.
6. Insert the spike by ¾ of its length into the outlet port. Do not touch the spike of the transfusion set and open the outlet port to avoid contamination.
7. After inserting the spike, squeeze the bag until the transfusion set filter is filled with blood.
8. Hook the bag from its holder. The system is ready for transfusion.

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