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KANSUK BLOOD AND BLOOD COMPONENTS PRESERVATION SOLUTIONS

CE1984

APPLICABLE CODES

20000006525	ACD-A ANTICOAGULANT SOLUTION 500 ML
-------------	-------------------------------------

DESCRIPTION

It is 500 ml ACD-A anticoagulant solution in transparent PVC bags, covered in a transparent package.

Each 1000 ml of ACD-A solution contains:

Glucose Monohydrate	24.5 g
Sodium Citrate Dihydrate	22.0 g
Citric Acid Anhydrous	7.3 g
Water for Injection q.s. ad	1000 ml


Figure-1: 20000006525

INTENDED USE

The ACD-A solution is used as an anticoagulant solution in apheresis devices during the collection of human blood per “Guideline for the Preparation, Use and Quality Assurance of Blood Components”.

PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

INDICATIONS FOR USE

This product is intended to be used for enabling producing blood components in automated blood collection devices (e.g. apheresis) by preventing coagulation for all population to aid any treatment. Blood components can be used for restoring blood components, improving oxygen delivery to tissues, restoring clotting factors and reducing or preventing bleeding according to the treatment decision of physician.

CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death. In addition, please see the contraindications of blood components produced by automated blood collection devices.

PRECAUTIONS

- The products should be stored in their original unopened / undamaged packaging at 0- +25 ° C, protected from direct sunlight.
- Use transfusion sets in compliance with EN ISO 1135-4.
- Follow validated blood component preparation procedures per local requirements.
- The tubes are suitable to be connected with a Sterile Connecting Device.
- It is intended for human blood and components.
- Not intended for direct intravenous infusion.
- Do not freeze.
- Comply with the apheresis protocol and the user manual of apheresis device.
- Patients or donors with impaired or abnormal citrate and/or calcium metabolism may show an increased risk of citrate sensitivity. The physician should evaluate the suitability of such patients or donors for apheresis and indicate how they should be monitored during the apheresis procedure.
- Single use only, after usage it should be thrown as a Medical Waste. **Caution:** Risk of infection when used more than once.
- Once the bag is opened, discard any unused solution.
- During shipment, up to two shipping boxes may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them, and they should be loaded in such a way to prevent them from falling over.

- For the explanation of symbols, see Symbols Description Table.

Caution: This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

GENERAL PREPARATIONS BEFORE USAGE

1. Open the outer package of the solution by pulling it from the corner and take the solution out of the package.
2. Check the bag for any visible damage. Do not use the defective sets.
3. Do not use if the solution is not clear. **Caution:** The solution may have a slightly yellow color. This is acceptable and does not affect safe use of the product.

Figure 2

PRACTICE

1. If the solution will be connected to a device, connect the solution according to the instructions of the device manufacturer.
2. Before transfusion, make sure that serological test results are acceptable and cross match test results with the recipient. If appropriate, perform the cross match test again at bed side.
3. Hold the top part of the T part with one hand and bag from the bottom part of T part with other hand, and twist of the upper part by 270° (Figure 2).
4. Keep holding the bag from T part, remove the cap of the transfusion spike with other hand.
5. Insert the transfusion spike of the apheresis set by ¾ of its length into the outlet port. Do not touch the spike of the transfusion set and open outlet port to prevent contamination (Figure 3).
6. After placing the spike, hang the bag to the apheresis device with the use of its hanger. ACD-A solution and apheresis set are ready to be used (Figure 4).

Figure 3

Figure 4

KANSUK LABORATUVARI SANAYİ VE TİCARET A.Ş.
Yassıören Mahallesi Fırat Sk. No:14/1 Arnavutköy-İstanbul/TURKEY
Tel: +90 212 592 15 76, Fax: +90 212 580 37 72
www.kansuk.com, e-mail: kansuk@kansuk.com

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