
I-TROMBOSİT HAVUZLAMA SETİNİN TANIMI

Ahtapot konfigürasyonu

REF20304046555 LABORATUVAR TROMBOSİT LÖKOREDÜKSİYON FİLTRE (IPP) / AHT

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DEFINITION OF THE I-PLATELET POOLING SET

Octopus configuration

REF 20304046555 LABORATORY
LEUKOREDUCTION FILTER for PLATELET,
(IPP) / OCT

One 700 ml pooling bag with 8 tubing leads for connection of buffy coats (BC) and one additional tubing lead for Platelet Additive Solution (PAS) dilution, an in-line Sepacell™ PLX-5 platelet leukodepletion filter connected to a 1,3 liter platelet storage bag with a 100 mL sampling pouch (Figure 1).


Figure-1

Train configuration

REF 20304146520 LABORATORY
LEUKOREDUCTION FILTER for PLATELET (IPP) /
TRAIN

A tubing lead with in line Sepacell™ PLX-5 platelet leukodepletion filter connected to a 1,3 liter platelet storage container with a 100 mL sampling pouch (Figure 2).


Figure-2

INTENDED USE

The I-Platelet Pooling Set is intended for the preparation of leukocyte depleted platelets obtained from whole blood derived pooled BC in additive solution. Validations were conducted using BC platelets prepared from pools of up to 8 BC.

- Up to 2 therapeutic doses of leukocyte depleted platelets can be obtained by pooling 7 or 8 BC and further treatment by a Pathogen Inactivation System (INTERCEPT™ Blood System).
- One therapeutic dose of leukocyte depleted platelets can be obtained by pooling 4 to 6 BC and further treatment by a Pathogen Inactivation System (INTERCEPT™ Blood System).
- One therapeutic dose of leukocyte depleted platelets can be obtained by pooling 4 to 6 BC without further Pathogen Inactivation treatment. Platelets can be stored at 22 + 2 °C under agitation up to 24 hours in the platelet storage container of the pooling set. The platelet concentrate must be transferred to another bag with infusion ports for transfusion.

PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

INDICATIONS FOR USE

This product is intended to be used for producing leukocyte-reduced platelets for all population to aid any treatment. Platelets can be used for restoring blood components and reducing and preventing bleeding according to the treatment decision of physician.

CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death. In addition, in case of a pathogen inactivation system is used, please see the user manual of the relevant system.

PRECAUTIONS

- A strict standardized and validated procedure for the preparation of venipuncture site should be in operation including skin disinfection and blood diversion of the first aliquots.
- Follow validated blood component preparation procedures per local requirements.
- In order to protect patients from an increased exposure (pools of > 6 BC) it is highly recommended that the leukocyte depleted platelets obtained with the pooling set are further processed immediately with a Pathogen Inactivation System (INTERCEPT™ Blood System) before the transfusion.
- In case of further processing with a Pathogen Inactivation System, strictly follow the instructions for use of the INTERCEPT™ Blood system.
- The tubes are suitable to be connected with a Sterile Connecting Device.
- It is intended for human blood and components.
- During shipment, up to two containers may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them.

and, they should be loaded in such a way to prevent them from falling over.

- Single use only, after usage it should be thrown as a Medical Waste. **Caution:** Risk of infection when used more than once.
- For the explanation of symbols, see Symbols Description Table.

CAUTION: This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. However, medical procedures should not be avoided as the benefits outweigh any potential health risks associated with exposure to DEHP. Please refer to the current literature to make an informed decision.

PROCESS STEPS

A. GENERAL PREPARATION BEFORE POOLING

1. For optimizing the platelet yield in the platelet concentrate, pre-selection of BC is recommended (donor platelet pre-count of at least $250 \times 10^9/L$).
2. Read manufacturer instruction for use for blood bags, centrifuge, blood component separator and procedures meeting local validation requirements to produce BC.
3. The recommended separation time for optimal filtration of the platelet concentrate is approximately 6 minutes. The blood component separator should be set up properly before use.

B. POOLING OF BUFFY COATS

Preparation before pooling

1. Select BC previously separated and stored.
2. Open the primary package and check the I- Platelet Pooling Set for any visible defects. Do not use defective sets.

Train method of pooling (indicative)

1. Sterile connect BC one to each other using the top and bottom tubings of the bags of the train.
2. Sterile connect the first BC of the train to a Platelet Additive Solution (PAS) bag.
3. Sterile connect the last (bottom) BC bag to the collection tube of the pooling set leaving 15-20 cm of tubing and close the clamp.
4. Place a re-usable clamp on the tube between the PAS bag and the first BC bag.
5. Hang the system from the PAS bag.
6. Let all BC flow by gravity into the last bag of the train.
7. To maximize the BC content recovery, rinse each BC bag sequentially using a portion of the PAS fluid. Let the remaining volume of PAS flow into the last bag of the train. Seal and detach the rest of the train from the last bag containing the pool.

Caution: Make sure that outlet ports are free from red blood cells (RBC) before centrifugation. RBC may pass into the platelet storage container or block the filter.

Octopus method of pooling (indicative)

1. Sterile connect BC (up to 8) to each dedicated tubing lead of the Octopus.
2. Close all the clamps.
3. Sterile connect a PAS bag to the last tubing lead of the Octopus.
4. Hang the set.
5. Open the clamps and let all the BC flow into the pooling bag.
6. Close the clamps.
7. To maximize the BC content recovery, rinse each BC bag sequentially using a portion of the PAS fluid. Let the remaining volume of PAS flow into the pooling bag. Seal and detach the Octopus from the pooling bag.

Caution: The pooling bag should be filled with a maximum of 650 ml of BC and PAS. Make sure that, the outlet ports are free from RBC before centrifugation. RBC may pass into the platelet storage container or block the filter.

C. CENTRIFUGATION

1. Pack the assembly of pooled buffy coats into centrifuge buckets according to blood center validated procedure.
2. Choose the appropriate platelet separation program. Follow user instructions of the centrifuge manufacturer.

Caution: Pack the sets into the centrifuge buckets carefully in order to avoid damages to the bags or poor separation results.

D. SEPARATION

1. After centrifugation, place the set on the blood component separator.
2. Hang the pooling bag onto the pins of the separator device.
3. Ensure the leukodepletion filter, tubings and platelet storage bag are placed on the device properly. Close the slide clamp of the sampling bag.
4. Break the breakaway cannula of the pooling bag. Start the device using the appropriate platelet separation program. Follow user instructions of the manufacturer.
5. The Press plate will press the collection bag and platelets will flow into the storage bag.
6. When RBC reach the tube downstream of the filter, stop the separation.
7. Seal and detach the platelet storage bag from the rest of the assembly.
8. To remove the air, gently fold the storage bag and squeeze the air into the sampling pouch. Once the air is removed, close the clamp and seal off the storage bag allowing to keep the 10-15 cm tubing for further processing with the pathogen inactivation system.

E. STORAGE

1. If further processing of the platelet concentrate is conducted with the INTERCEPT™ Blood System, platelets must be stored at 22 ± 2 °C under continuous agitation. Strictly follow the instructions for use of the INTERCEPT™ Blood System for further processing
2. If no further processing with a Pathogen Inactivation System will be performed, the leukocyte depleted platelets obtained with the pooling set should be used within 24 hours from the platelet pool production.

KANSUK LABORATUVARI SANAYİ VE TİCARET A.Ş.
Yassören Mahallesi Fırat Sk. No:14/1 Arnavutköy-İstanbul/Turkey
Tel: +90 212 592 15 76, Fax: +90 212 580 37 72
www.kansuk.com, E-mail: kansuk@kansuk.com

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