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20004046520

LABORATUVAR TROMBOSİT LÖKOREDÜKSİYON FİLTRE PRP

TANIM

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APPLICABLE CODES

20004046520	LABORATORY LEUKOREDUCTION FILTER for PLATELET, PRP
-------------	--

DESCRIPTION

Laboratory type platelet filter is composed of a 1300 ml platelet storage bag, a platelet filter connected to the top of the tube which is at the end of the storage bag and a collection tube connected to the filter outlet, an aggregate filter on the tube and a spike. It consists of a bypass line for removing air inside the bag (Figure 1 and 2). As optional, in some systems there is a pH sensor to be connected to the pH SAFE reader device to read pH in the storage bag.


Figure 1


Figure 2

INTENDED USE

It is used for obtaining one therapeutic dose of leukocyte depleted platelet suspension by pooling 4-6 units of "Random Platelet Suspension".

PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

INDICATIONS FOR USE

This product is intended to be used for producing leukocyte-reduced platelets for all population to aid any treatment. Platelets can be used for restoring blood components and reducing and preventing bleeding according to the treatment decision of physician.

CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

PRECAUTIONS

- Use transfusion sets in compliance with EN ISO 1135-4.
- Follow validated blood component preparation procedures per local requirements.
- The tubes are suitable to be connected with a Sterile Connecting Device.
- It is intended for human blood and components.
- During shipment, up to two containers may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them, and, they should be loaded in such a way to prevent them falling over.
- Single use only, after usage it should be thrown as a Medical Waste. **Caution:** Risk of infection when used more than once.
- For the explanation of symbols, see Symbols Description Table.

Caution: This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

PROCESS STEPS

A. GENERAL PREPARATION BEFORE POOLING

1. Open the primary package, remove the laboratory type platelet filter and check it for any visible defects. Do not use defective sets.
Caution: The bag system may give off a different smell because of the gamma sterilization.
2. Prepare 4-6 units of the random platelets from the same group.
3. Set up the filter hanger to enable the connected bags to pool by filtering.

B. FILTRATION

1. Weld one of the selected random platelets to the outlet edge of the filter in the collection bag using the sterile connection device. Please read SCD device manual for its usage. You can enter the random platelet bags using a spike. In this case, this process must be done under laminar flow.
2. Hang the platelet bag to the filter hanger.
3. Collect the platelets filtered into the pooling bag.
4. For the remaining random platelets, follow the same procedure starting from the item 2.
Caution: When the second and subsequent bags are being connected, you must be sure that the front aggregate filter is filled with at least 2 cm of platelet suspension to prevent the filter from taking air in.
5. After transferring all the selected platelets to the collection bag, close the temporal clamp under the filter and release the trapped air in the collection bag to the empty platelet bag above or to the sampling bag using the bypass system which is on the same line.
6. Close the collection tube by the tube sealing device at a suitable point, 15-20 cm above after the tube outlet and remove from the system.
7. One therapeutic dose of leukocyte depleted platelet suspension is obtained.
8. Please see the storage section for storage conditions.

C. STORAGE

1. Platelets can be stored at 22 °C under agitation up to 5 days in a platelet bag made of TOTM.
2. In containers suitable for storage up to 7 days, platelets can be stored up to 7 days provided that regular invasive controls of bacterial contamination are performed.
3. In containers equipped with pH sensor, by performing regular non-invasive controls of pH, the status of platelets can be determined. **Caution:** Please read BCSI pH1000 user manual instructions.

D. TRANSFUSION

1. Inspect the container and/or blood component for any defects.
2. Before transfusion, make sure that serological test results are acceptable and cross match test results with the recipient. If appropriate, apply cross match again at bedside.
3. Platelets should be inspected before transfusion to ensure that no platelet aggregates are visible.
4. Hold the top part of the T part with one hand and bag from the bottom part of T part with other hand, and twist of the upper part by 270°.
5. Keep holding the bag from T part, remove the cap of the transfusion spike with other hand.
6. Insert the spike by ¾ of its length into outlet port. Do not touch the spike of the transfusion set and open outlet port to prevent contamination.
7. After inserting the spike, squeeze the bag to the level that transfusion set filter is filled with blood.
8. Hang the bag onto the stand. The unit is ready for transfusion.

KANSUK LABORATUARI SANAYİ VE TİCARET A.Ş.
Yassıören Mahallesi Fırat Sk. No:14/1 Arnavutköy-İstanbul/Turkey
Tel: +90 212 592 15 76, Fax: +90 212 580 37 72
www.kansuk.com, E-mail: kansuk@kansuk.com

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