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
APPLICABLE CODES

68135	LABORATORY LEUKOREDUCTION FILTERS for RBC
20014068130	LABORATORY LEUKOREDUCTION FILTER for RBC TYPE I
68137	LABORATORY LEUKOREDUCTION FILTER for RBC TYPE II

DESCRIPTION

They are disposable, sterile transfusion sets which contain a transfer spike, di-(2-ethylhexyl)-phthalate (DEHP) plasticized polyvinyl chloride (PVC) tubing, a blood bag, a leukocyte filter and drip chamber (Figure 1 and Figure 2). They are sterilized by Gamma irradiation.


Figure 1 - 20014068130

Figure 2 - 68137


Figure 3 - 68135

INTENDED USE

For preparation or storage of leukodepleted red blood cell (RBC) suspensions per "Guide to the Preparation, Use and Quality Assurance of Blood Components".

PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

INDICATIONS FOR USE

This product is intended to be used for leukocyte reduction of red blood cells for all population to aid any treatment. Red blood cells can be used for restoring blood components and improving oxygen delivery to tissues according to the treatment decision of physician.

CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

PRECAUTIONS

- The tubes are suitable to be connected with a Sterile Connecting Device.
- It is intended for human blood and components.
- During shipment, up to two containers may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them, and, they should be loaded in such a way to prevent them from falling over.
- Single use only, after usage it should be thrown as a Medical Waste. **Caution:** Risk of infection when used more than once.

- For the explanation of symbols, see Symbols Description Table.

Caution: This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. However, medical procedures should not be avoided as the benefits outweigh any potential health risks associated with exposure to DEHP. Please refer to the current literature to make an informed decision.

PROCESS STEPS

A. GENERAL PREPARATIONS BEFORE USE

1. Take the filter out from primary packaging. Filter is ready to be connected to the blood bag (Figure 1 and Figure 2).
2. Please check the filter set for any visible defects. Do not use defective sets.
3. Filter should be connected to the bag on a proper bench.

Caution: If you use the system showed in Figure 1 and Figure 2, make sure the clamp under the drip chamber is closed (Figure 1 and Figure 2 no.8).

B. FILTRATION

1. Connect the filter system to the bag to be filtered.
 - a. If connection will be performed by the sterile connector device (SCD), the tube of the RC unit and tube of the filter set have to be connected (Figure 1 and Figure 2, No: 3e and Figure 3). The tubing of the set and the tubing of the RC unit should be inserted in the SCD with tubings looking out of the SCD device about 4-5 cm.

Caution: See user manual of the SCD.

Figure 4

- b. If the connection will be performed by spiking (Figure 4);
 - i. Open one of the outlet ports of the bag. **Caution:** See blood bag user manual.
 - ii. Open the safety cap of the spike of the filter (Figure 1, No: 5) and connect it to the outlet port of the blood bag firmly.

Caution: Make sure that the air vent cap of the spike is closed.

Figure 5

2. Mix the red blood cell suspension slowly and hang the blood bag (Figure 6).

Figure 6

3. Open the clamp under the drip chamber (Figure 1 and Figure 2, No: 8, and Figure 7).

Figure 7

4. Firstly, ¼ of drip chamber is automatically filled with blood after opening the clamp and then the blood flows through the leukocyte filter (Figure 1 no.12 and Figure 2 no.12).

Caution: Do not touch the drip chamber and observe the filter is primed properly for 1 min (Figure 8)

Figure 8

5. Firstly, blood flows through the front side of the leukocyte filter (Figure 8). After that, the back side of the filter is filled with blood (Figure 9) and leukocyte-depleted blood flows through into the storage bag.

Figure 9


Figure 10

6. After all of the blood is filtered, close the clamp under the filter (Figure 1 no: 8, Figure 2 no:8 and Figure 11).

Figure 11

7. Squeeze the transfer bag (Figure 12) and push the air in the bag out through the check valve (Figure 1 No:3b and Figure 2 No:3b) and transfer the air into the main bag. Then re-open the clamp and wait for the remaining blood flows through the filter and tubing until filter turns pink (Figure 13).

Caution: Preferred air removal process is done while sitting as the angle of air removal is much better and easier as compared to a standing position.

Figure 12

Figure 13

8. The median filtration time is around 15-20 minutes depending on hematocrit and temperature.
9. When all the blood is filtered completely, remove the filtered RCs from the system by sealing its tubing with a sealer, it will then be ready for storage or usage.

C. STORAGE

- Red blood cells should be stored under 1-6 °C at blood bank refrigerators.

Caution: While removing blood bags from the blood bank refrigerator in case of need, it is very important to handle them carefully so that they are not exposed to hard and sharp objects that could damage them and not rubbed to the floor and dropped or thrown from a high distance while placing them on the table.

Caution: After removing the blood bags from the blood bank refrigerator, they should be kept at room temperature for 4-10 minutes to return to the normal temperature since the PVC structure of the blood bag has a fragile nature at colder temperatures.

D. TRANSFUSION

1. Use transfusion sets in compliance with EN ISO 1135-4.
2. Inspect the container and/or blood component for any defects.
3. Before transfusing the blood, make sure that serological test results are acceptable, cross match test results match with the recipient. If appropriate, apply cross match again at bedside.

*** Please read steps 4-5-6-7-8 carefully if your bag type is as shown in the picture.

4. Hold the top part of the T part with one hand and bag from the bottom part of T part with other hand, and twist of the upper part by 270°
5. Keep holding the bag from T part, remove the cap of the transfusion spike with other hand.
6. Insert the spike by $\frac{3}{4}$ of its length into outlet port. Do not touch the spike of the transfusion set and open outlet port to prevent contamination.
7. After inserting the spike, squeeze the bag to the level that transfusion set filter is filled with blood.
8. Hook the bag from its holder. The system is ready for transfusion.

Figure 14

*** Please read steps 9-10-11-12-13 carefully if your bag type is as shown in the picture.

9. While holding the bottom of the auricle (OUTLET PORT) with one hand, grab the top of the auricle with the other hand and tear the auricle from the upper section.
10. Remove the cap of the transfusion pin with the other hand while continuing to hold the bag by the auricle.
11. Insert $\frac{3}{4}$ of the tip of the transfusion set into the outlet port. To avoid contamination, do not touch the spike of the transfusion set and the exposed outlet port.
12. After inserting the pin, tighten the bag until the filter of the transfusion set is filled with blood.
13. Hook the bag from its holder. The system is ready for transfusion.

Figure 15

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