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
APPLICABLE CODES

20011268115	BED-SIDE LEUKOREDUCTION FILTER for RBC TYPE I
20014368115	BED-SIDE LEUKOREDUCTION FILTER for RBC TYPE III
20014068115	BED-SIDE LEUKOREDUCTION FILTER for RBC
20013068115	BED-SIDE LEUKOREDUCTION FILTER for RBC TYPE II

DESCRIPTION

They are disposable, sterile transfusion sets which contain a transfer spike, di-(2-ethylhexyl)-phthalate (DEHP) plasticized poly(vinyl chloride) (PVC) tubing, a leukocyte filter, drip chamber, a roller clamp and luer connection. (Figure 1a, 1b, 1c, 1d). They are sterilized by Gamma irradiation.


Figure 1a

20011268115


Figure 1b

20014368115


Figure 1c

20014068115


Figure 1d

20013068115

INTENDED USE

For leukodepletion of red blood cell suspensions and transfusion at bedside per "Guide to the Preparation, Use and Quality Assurance of Blood Components".

PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by nurses in the clinics.

INDICATIONS FOR USE

This product is intended to be used for leukocyte reduction and transfusion of red blood cells for all population to aid any treatment. Red blood cells can be used for restoring blood components and improving oxygen delivery to tissues according to the treatment decision.

CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

PRECAUTIONS

- The tubes are suitable to be connected with a Sterile Connecting Device.
- It is intended for human blood and components.
- During shipment, up to two containers may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them, and, they should be loaded in such a way to prevent them from falling over.
- Single use only, after usage it should be thrown as a Medical Waste. **Caution:** Risk of infection when used more than once.
- For the explanation of symbols, see Symbols Description Table.

CAUTION: This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. However, medical procedures should not be avoided as the benefits outweigh any potential health risks associated with exposure to DEHP. Please refer to the current literature to make an informed decision.

PROCESS STEPS

A. GENERAL PREPARATION BEFORE USE

1. Take the filter out from primary packaging. Filter is ready to be connected to the blood bag.
2. Please check the filter set for any visible defects. Do not use defective sets.
3. Close the clamp under the spike (Figure 2).
4. Filter should be connected to the bag on a proper bench.

Figure 2

B. FILTRATION AND TRANSFUSION

1. Connect the filter system to the bag to be filtered.
 - a. If the connection will be performed using sterile connecting device, the tube of the RC unit and the tube of the filter set have to be connected (under the spike) using sterile connecting device (Figure 3). **Caution:** See user instructions for sterile connecting device.

Figure 3

- b. If the connection will be performed by spiking;
 - i. Open one of the outlet ports of the bag. **Caution:** See instructions for use blood bag.
 - ii. Open the safety cover of the spike of the filter and connect to the outlet port of the blood bag without letting the blood to leak. (Figure 4) **Caution:** Make sure that the air vent cap of the spike is closed.

Figure 4

2. Hang the blood bag (Figure 5).

Figure 5

3. Open the clamp under the spike.
4. Blood flows into the filter (Figure 6).

Figure 6

5. After the filter is filled with blood, invert the drip chamber and fill it half way with blood. Then, put the drip chamber into its normal position. In this way, the air present in the system is removed. (Figure 7)


Figure 7

6. Wait until the blood reaches to Luer Lock and then stop the blood flow by closing the roller (Figure 8).

Figure 8

7. Open the luer lock cap, connect the needle (Figure 9) and transfuse to patient.

Caution: Inspect the bag and blood for any defects.

Caution: Before transfusing the blood, make sure that serological test results are acceptable, cross match test results match with the recipient. If appropriate, apply cross match again at bedside.

Figure 9

8. Adjust the blood flow by roller. Suggested transfusion time and speed: 19-20 mmHg drop/ min.

9. To let the blood in the drip chamber with pre-filter flow after the blood in the bag is removed or to let the blood in the bag flow after the 10-30 ml of blood remained in the bag;

- a. For the sets shown in Figure 1a, 1b and 1c open the cap of the vented spike (Figure 10).

Figure 10

- b. For the set shown in Figure 1a and 1b, open the air valve (No:12 and 4) under the filter after the blood in the bag and pre-filter drip chamber is removed (Figure 11).

Figure 11

10. After all process, transfusion of one unit red blood cell suspension to patient is already completed.

KANSUK LABORATUVARI SANAYİ VE TİCARET A.Ş.
Yassören Mahallesi Fırat Sk. No:14/1 Arnavutköy-İstanbul/Turkey
Tel: +90 212 592 15 76, Fax: +90 212 580 37 72
www.kansuk.com, E-mail: kansuk@kansuk.com

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