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
### APPLICABLE CODES

20004068145	BED-SIDE LEUKOREDUCTION FILTER for PLATELET
-------------	---

### DESCRIPTION

It is a disposable, sterile transfusion set which contain a transfer spike, di-(2-ethylhexyl)-phthalate (DEHP) plasticized poly(vinyl chloride) (PVC) tubing, a leukocyte filter, drip chamber, a roller clamp and luer connection.

### INTENDED USE

For leukodepletion and transfusion of platelets (PLTs) obtained by platelet rich plasma at bedside per “Guide to the Preparation, Use and Quality Assurance of Blood Components”.

### PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

### INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by nurses in the clinics.

### INDICATIONS FOR USE

This product is intended to be used for leukocyte reduction and transfusion of platelets for all population to aid any treatment. Platelets can be used for restoring blood components and reducing and preventing bleeding according to the treatment decision.

### CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

### PRECAUTIONS

- The tubes are suitable to be connected with a Sterile Connecting Device.
- It is intended for human blood and components.
- During shipment, up to two containers may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them, and, they should be loaded in such a way to prevent them from falling over.
- Single use only, after usage it should be thrown as a Medical Waste. **Caution:** Risk of infection when used more than once.
- For the explanation of symbols, see Symbols Description Table.

**CAUTION:** This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

## **PROCESS STEPS**

### **A. GENERAL PREPARATION BEFORE USE**

1. Take the filter out from primary packaging. Filter is ready to be connected to the platelet bag (Figure 1).
2. Please check the filter set for any visible defects. Do not use defective sets.
3. Select and prepare 4-6 units of random PLTs to be transfused to the patient.
4. Filter should be connected to the platelet bag on a proper bench.

### **B. FILTRATION AND TRANSFUSION**

1. Close the clamp under the spike (Figure 3).


Figure 3

2. Open the safety cap of the filter spike (Figure 1, No: 9) and connect it to the outlet port of the platelet bag without letting the blood to leak (Figure 4).

Figure 4

3. Hang the random PLT bag connected to filter at bed side. (Figure 5)

Figure 5

4. Open the clamp under the spike.
5. Once the filter is filled, PLTs begin to fill up the drip chamber (Figure 6)

Figure 6

6. Wait until the PLTs reach to the luer lock and then stop the flow by closing the roller clamp (Figure 7).

Figure 7

7. Open the luer lock cap, connect the luer lock to the scalp etc. that is put into the vascular access of the patient before (Figure 8).

**Caution:** Inspect the bag and blood for any defects.

**Caution:** Before transfusing the blood, make sure that serological test results are acceptable, cross match test results match with the recipient. If appropriate, apply cross match again at bedside.

Figure 8

8. Open the roller clamp and adjust the blood flow. Recommended transfusion rate is 260-320 mL/h.
9. When 10-20 ml of PLTs has been remained in the PLT bag, close the clamp under the spike and pull out the spike of filter from the PLT bag and connect it to a new PLT bag.
10. Apply this process for all random platelets to be connected to the patient in the same way.

**NOTE:** When 10-20 ml of PLTs have remained in the random platelet bag, close the clamp and connect the spike of the filter to the new platelet bag without contacting with air. Otherwise, when all the platelets in the bag are over, air may pass into the filter and cause obstruction in the filter and the platelets inside the next platelet bag may not pass through the filter. Give all of the platelet to be last transfused the patient.

**KANSUK LABORATUVARI SANAYİ VE TİCARET A.Ş.**  
Yassıören Mahallesi Fırat Sk. No:14/1 Arnavutköy-İstanbul/Turkey  
Tel: +90 212 592 15 76, Fax: +90 212 580 37 72  
[www.kansuk.com](http://www.kansuk.com), E-mail: [kansuk@kansuk.com](mailto:kansuk@kansuk.com)

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