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
20000006217 | KAN TORBASI TEKLİ CPDA-1 450 ML

TANIM

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
DESCRIPTION

Closed blood bag systems made of PVC, consisting of a 16-g needle, a blood bag and including CPDA-1 as the anticoagulant solution.

INTENDED USE

For collection and storage of whole blood per "Guide to The Preparation, Use and Quality Assurance of Blood Components."

PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

INDICATIONS FOR USE

This product is intended to be used for collection of whole blood for all population to aid any treatment. Whole blood and blood components derived from whole blood can be used for restoring blood components, improving oxygen delivery to tissues, restoring clotting factors and reducing or preventing bleeding according to the treatment decision of physician.

CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

PRECAUTIONS

- Do not use if needle protector is opened.
- Do not add medication to blood.
- Cross match the contents of the bag with recipient before transfusion.
- Use transfusion sets in compliance with EN ISO 1135-4.
- Follow validated blood component preparation procedures per local requirements.
- The tubes are suitable to be connected with a Sterile Connecting Device.
- It is intended for human blood and components.
- During shipment, up to two shipping boxes may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them, and, they should be loaded in such a way to prevent them from falling over.
- Single use only, after usage it should be thrown as a Medical Waste. **Caution:** Risk of infection when used more than once.
- For the explanation of symbols, see Symbols Description Table.


Figure 1

Caution: This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

GENERAL PREPARATIONS BEFORE USAGE

- Open the outer wrap at the tear notches and remove one blood bag system, check for any visible defects of tubes, needle, anticoagulant and additive solutions in terms of appropriate volume, discoloration or particulate contaminants.

Caution: Do not use defective sets, make sure the solutions are clear.

- Apply blood pressure cuff or tourniquet on donor's arm, identify venipuncture site and release the cuff or tourniquet.
- Disinfect the venipuncture site. Cover the area with sterilized gauze strip and do not touch until the time of venipuncture.

Caution: Beware of iodine sensitive donors!

- Signs on the needle protector and hub should be the in the same line (Figure 1). Apply some pressure to determine any leakages.
- Place the collection bag below the level of donor's arm.

Caution: While opening the blood bag system, blood taking tube should be held from needle protector not from the needle safety cover. Pay attention not to lock the needle safety cover before collection of blood by mistake.

BLOOD COLLECTION (Indicative)

- Reapply tourniquet or inflate blood pressure cuff until previously selected vein is prominent again.
- Hold the needle from marked side. To remove the needle cap, twist to break the seal between the needle and its cap and remove it by sliding it down the axis of the needle. **Caution:** Be sure not to touch the cannula!
- Apply the venipuncture immediately. Venipuncture is very important for a non-coagulated blood.
- Strap the needle to the donor's arm.
- Mix blood and anticoagulant gently and periodically (approximately in every 45 seconds) during collection. If a mixer or a collection system is used, please refer to their operator's manual.

6. Make sure there is continuous blood flow. If there is continuous, adequate blood flow and constant agitation, blood collection is expected to be completed under 12 minutes. **Caution:** If it is above 12 minutes, collected blood may not be suitable for Platelets, Fresh Frozen Plasma or Cryoprecipitate AHF.
7. Collect the quantity of blood within the limits indicated on the bag label. Monitor the blood being collected.
8. After blood collection is completed, stop the blood mixing device and loosen the tourniquet or deflate the pressure cuff.
Caution: To prevent accidental injuries or contamination of the personnel, insert the needle into the locked system or throw it to the medical waste.
9. Close the blood taking tube via tube sealing device and remove the needle from the system. **Caution:** Throw locked needle kit to the medical waste.
10. Strip donor tubing as much as possible through the bag and start sealing. Work quickly to prevent the blood from clotting in the tubing. Invert bag several times to mix thoroughly, then allow tubing to refill with anticoagulated blood from the bag. Repeat this procedure twice. (Figure 2).
11. Check the blood bag again for any faults.
12. Make sure all donor information required is filled out on the label on the blood bag.
13. After collection, whole blood should be placed immediately at 1-6 °C in a blood bank refrigerator which is designed for this purpose to be stored up to 35 days.

Figure 2

TRANSFUSION

1. Use transfusion sets in compliance with EN ISO 1135-4.
2. Inspect the container and/or blood component for any defects.
3. Before transfusing the blood, make sure that serological test results are acceptable, cross match test results match with the recipient. If appropriate, apply cross match again at bedside.

*** Please read steps 4-5-6-7-8 carefully if your bag type is as shown in the picture.

4. Hold the top part of the T part with one hand and bag from the bottom part of T part with other hand, and twist of the upper part by 270°
5. Keep holding the bag from T part, remove the cap of the transfusion spike with other hand.
6. Insert the spike by ¾ of its length into outlet port. Do not touch the spike of the transfusion set and open outlet port to prevent contamination.
7. After inserting the spike, squeeze the bag to the level that transfusion set filter is filled with blood.
8. Hook the bag from its holder. The system is ready for transfusion.

Figure 4

*** Please read steps 9-10-11-12-13 carefully if your bag type is as shown in the picture.

9. While holding the bottom of the auricle (OUTLET PORT) with one hand, grab the top of the auricle with the other hand and tear the auricle from the upper section.
10. Remove the cap of the transfusion pin with the other hand while continuing to hold the bag by the auricle.
11. Insert ¾ of the tip of the transfusion set into the outlet port. To avoid contamination, do not touch the spike of the transfusion set and the exposed outlet port.
12. After inserting the pin, tighten the bag until the filter of the transfusion set is filled with blood.
13. Hook the bag from its holder. The system is ready for transfusion.

Figure 3

Figure 5

KANSUK LABORATUVARI SANAYİ VE TİCARET A.Ş.
Yassiören Mahallesi Fırat Sk. No:14/1 Arnavutköy-İstanbul/Turkey
Tel: +90 212 592 15 76, Fax: +90 212 580 37 72
www.kansuk.com, e-mail: kansuk@kansuk.com

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