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2210006732

KAN TORBASI T&amp;B ERİTROSİT LÖKOSİT FİLTRELİ DÖRTLÜ CPD/SAGM 450 ML

## TANIM

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**APPLICABLE CODES**

22100006732	BLOOD BAG T&B WITH RED BLOOD CELL LEUKOCYTE FILTER QUADRUPLE CPD/SAGM 450 ML
-------------	--

**DESCRIPTION**

Closed blood bag system made of PVC consisting of a erythrocyte leucocyte sampling system, 16 g needle, transfer bags and including CPD as anticoagulant solution and SAG-M as additive solution.

**INTENDED USE**

Used for leukocyte-reduced erythrocyte suspension and plasma and platelet collection, separation and storage.

**PATIENT POPULATION**

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

**INTENDED USER AND INTENDED USE ENVIRONMENT**

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

**INDICATIONS FOR USE**

This product is intended to be used for collection of whole blood and preparing leukocyte-reduced blood components from whole blood for all population to aid any treatment. Blood components can be used for restoring blood components, improving oxygen delivery to tissues, restoring clotting factors and reducing or preventing bleeding according to the treatment decision of physician.

**CONTRAINDICATIONS**

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

**GENERAL PREPARATIONS BEFORE USAGE**

8. Open the outer wrap of the blood bag through the tearing notch and remove a blood bag system.
9. Attach the sphygmomanometer or tourniquet on donor's arm, identify venipuncture site and release the tourniquet or sphygmomanometer.
10. Disinfect the venipuncture site. Cover the area with sterilized gauze strip and do not touch until the time of venipuncture.

**CAUTION:** Beware of iodine sensitive donors!

11. Before entering the vein, review the blood bag system, tubings and needle for visible errors and color change.
12. Signs on the needle protector and hub should be in the same line. Apply some pressure to determine any leakages.
13. Place the collection bag below the level of donor's arm.
14. Check the anticoagulant solution and additional solutions for the correct volume, color and particle content.

**CAUTION:** While opening the blood bag system, blood taking tube should be held from needle protector not from the needle safety cover. Pay attention not to lock the needle safety cover before collection of blood by mistake.

**INSTRUCTION FOR USE**

**BLOOD COLLECTION:**

**IF BLOOD BAG SYSTEMS WITH SAMPLING ARE USED, SEE POINTS 1-4; IF THE SYSTEM DOES NOT HAVE SAMPLING, IF BLOOD BAG SYSTEMS WITH SAMPLING ARE USED, SEE POINTS 1-4; IF THE SYSTEM DOES NOT HAVE SAMPLING, CONTINUE FROM POINT 5.**

1. Before removing the needle's cap, close the clamp on the tubing to collection bag (for systems in Figure 1 and 2). If the system's donor line has breakaway valves, not necessary to close the clamp off.
2. Break open the breakaway valve on the tubing to the sampling holder / bag (for systems in Figure 1 and 2). **NOTE:** There are not any breakaway valves on the tubing to the sampling tubing holder / bag in the systems with a blood collection line with breakaway valves.
3. Depending on the sampling system used, follow one of the following processes:
  - c. **If a system with tube holder is used;**
    - i. Insert the sampling tubes into the tube holder and collect the required blood sample (Figure 1, No: 3).
    - ii. **CAUTION:** First, take the samples into bigger tubes.
    - iii. Open the clamp on the tubing to the bag and start blood collection (Figure 1, No: 1).
  - d. **If a system with sampling bags is used;**
    - i. Take the blood to the sampling bag (Figure 2, No: 5).
    - ii. Close the clamp on the tubing to the sampling bag (Figure 2, No: 2).
    - iii. Open the clamp on the tubing to the bag and start blood collection (Figure 2, No: 1).
    - iv. Place the sampling bag in the palm tube holder facing down. Insert the sampling tubes into the tube holder connected to the sampling tube and collect the required blood sample (Figure 2, No: 4). **CAUTION:** Blood sample in the sampling bag should immediately transferred to the tubes; otherwise, the blood may get coagulated.
4. After the sampling process is finished, break the breakaway valve on the tubing to the sampling tubing bag in the systems with a blood collection line with breakaway valves (Figure 3 and 4) and open the clamp in others (Figure 1 and 2).
5. Hang the tube with the tube holder onto donor's wrist in a way not to compress the needle. Continue from point 9.
6. If the system used does not have sampling, attach a temporary clamp or knot the blood collection tubing.
7. Until the previously selected vein is visible again, inflate the sphygmomanometer again or apply the tourniquet.
8. Hold the needle from marked side. To remove the needle casing, turn the cap in a circular motion so that the adhesion zone between the needle and the casing is broken and scrape down along the needle axis. **CAUTION:** Be sure not to touch the needle!
9. Apply the venipuncture with the needle immediately. It is very important to apply the venipuncture cleanly and smooth for a complete blood collection without clots.
10. Strap the needle to the donor's arm.
11. Mix blood and anticoagulant gently and periodically (approximately in every 45 seconds) during collection, If a mixer or a collection system is used, please refer to their operator's manual.
12. Make sure there is continuous blood flow. If there is continuous, adequate blood flow and constant agitation, blood collection is expected to be completed within 12 minutes.

**CAUTION:** If it is above 12 minutes, collected blood may not be suitable for Platelets, Fresh Frozen Plasma or Cryoprecipitate AHF separation.

13. Collect the quantity of blood within the limits indicated on the bag label. Monitor the blood being collected.
14. When the blood collection process is complete, the blood collection shaker takes the blood at the set weight or milliliter; then stop the blood collection by attaching a clamp to the blood collection tubing.
15. After blood collection is completed stop the blood mixing device and loosen the tourniquet or deflate the pressure cuff.
16. The clamp of the blood collection shaker is turned on or lower the pressure of the sphygmomanometer.
17. When the blood collection process is completed, to remove the needle, hold it with the thumb and index fingers of the other hand while holding the sterile gauze in one hand, and remove the needle completely without applying pressure. Let the donor gently press the sterile gauze to the venipuncture point.
18. Holding the needle with one hand, wait for the blood from the blood collection tubing to flow into the bag so that the needle tip is upwards.
19. After this process, close the clamp on the tubing to the bag. In systems without sampling, knot the tubing. **CAUTION:** If the clamp does not close, there may be a risk of blood gushing out of the needle protection cover while the needle is locked.
20. For the blood bag systems with needle safety cover, while holding the needle with one hand hold the needle safety cover with the other one and then leave the needle, pull the needle to the safety cover by pulling the tubing under the safety cover and be sure that it is locked (Figure 5, No: 6). **CAUTION:** To prevent accidental injuries or contamination of the personnel in case of systems without needle safety cover, insert the needle into the needle destruction machine or throw it to the bio-danger box.
21. Close the blood taking tubing under the Y part by tubing sealing device and remove the needle from the system. **CAUTION:** Throw locked needle system to the medical waste.
22. Scrape the blood from the blood collection tubing completely from the place of origin to the bag. Work quickly to prevent the blood from clotting in the tubing. Stir the bag by inverting it over a few times, then allow tubing to refill with anticoagulated blood. Repeat this procedure twice (Figure 6).
23. For cross-comparison purposes, segment the tubing connected to the collection bag with the tubing shut-off machine so that the segment numbers are clearly and fully readable. The segment closest to the collection bag is always attached to the bag for identification purposes.
24. Check the blood bag again for any faults.
25. Make sure that all donor information required is filled out on the label on the blood bag.
26. After complete blood collection, immediately place it in a blood bank cabinet designed for this process at 1-6°C. If platelets are to be separated, do not cool the blood, but leave it at an angle of 45° until the temperature reaches 20-24°C (min one hour). For gaining fresh frozen plasma and platelets, complete all the steps within 8 hours after collection.


Figure 1 Figure 2 Figure 3 Figure 4


Figure 5

Figure 6


Figure 7

## USE OF BAGS:

### TOP&BOTTOM 3- AND 4-PIECE BAG SYSTEMS (Used for leukocyte-reduced erythrocyte, plasma, platelet separation according to buffy coat (BC) method)

1. Load the bag system into the centrifugal rods as explained in point 1 so that the breakaway valve of the main bag is looking up.
- CAUTION:** While loading the godes; IF HAS A HARD FILTER; Place the filter in the filter compartment by placing it in the filter compartment. IF IT HAS A SOFT FILTER; Place the filter neatly between the transfer bags and place the whole system in the centrifuge gode. Proper positioning of the bags will prevent burst of the bags during centrifugation.
2. If the blood taken is not separated into components but is to be used as whole blood, break the breakaway valve on the main bag and transfer the whole blood to the transfer bag with a capacity of 500 ml in the separator.
3. Centrifuge the entire system at 3500 g for 15 minutes under 22°C. **CAUTION:** Centrifuge settings may vary according to the centrifuge brand used.
4. After the centrifuge, place the main bag in the automatic separator device. If quad system is used, close the tubing of the platelet bag with a clamp.
5. Use the automatic separator according to the instructions for use.
6. First, break the breakaway valve of the main bag, then the breakaway valves of the other transfer bag (Figure 7).
7. The automatic separator automatically transfers the plasma to the plasma bag and the erythrocytes to the transfer bag, with the Buffy Coat remaining in the main bag.
8. When the process is completed, close the transfer bags with plasma and erythrocyte suspension with a tubing closure device and disconnect them from the system.
9. Store the plasma in the freezer. If the separation takes place within 6 hours after collection, it can be stored at -18°C/-25°C for 3 months, -30°C/-40°C for 6 months and at lower temperatures for 12 months.
10. Hang the SAG-M bag on the filter hanger at a height of 150 cm from the ground (Figure 8).
11. Break the breakaway valve at the opening of the bag by turning left and right (Figure 7).
12. The SAG-M solution passes through the filter by the force of gravity.
13. When the SAG-M transition process is finished, pressurize it with a temporary clamp under the leukocyte filter on the hose.
- CAUTION:** Make sure that the clamp is made from the point closest to the erythrocyte bag.
14. Gently stir the erythrocyte bag with SAG-M added on it before hanging on the hanger. **CAUTION:** Performing the mixing process shortens the filtration time.
15. Flip the entire system upside down, hang the erythrocyte bag on the top and leave the SAG-M bag below (Figure 9).
16. Turn on the temporary clamp open. With the force of gravity, erythrocytes pass through the filter and filter the leukocytes and collect them in the SAG-M bag.
- CAUTION:** In order for the filter not to make air and to increase the filtration rate, make sure that the tubing is flat after hanging the bags, pull the clamp at once.
17. The filtration time is 25 minutes on average at controlled room temperature.
18. Scrape the blood from the blood collection tubing completely from the place of origin to the bag and make sure the tubing is filled with anticoagulated blood. Repeat this procedure twice (Figure 6).
19. RBCs should be placed into blood bank refrigerators designed for this process under 1-6°C.
20. If a quadruple bag is used, close the middle of the main bag with the buffy-coat with a platelet clamp holder, mix the bag well and leave for 1 hour.
21. Centrifuge the main bag and 5 inscribed transfer bags at 500 g for 6 min. at 20°C. **CAUTION:** Centrifuge settings may vary according to the centrifuge brand used.
22. Transfer platelets to the platelet bag in the automatic separator. **CAUTION:** After centrifuge, platelets are immediately separated.
23. Platelets can be stored in the agitator with incubator for up to 5 days inside the transfer bag made of special PVC (TOTM).

Figure 8


Figure 9

## BLOOD DONATION

1. Inspect the container and blood for any faults.
2. Before applying the blood, check whether the serological tests applied are appropriate and whether the results of the cross-comparison test are consistent with the blood given. If appropriate, apply cross match again at bedside.
3. Use a transfusion kit in accordance with EN ISO 1135-4.
4. Hold the top part of the T part with one hand and bag from the bottom part of T part in the bag with other hand, and twist of the upper part by 270° (Figure 10).
5. Remove the cap of the transfusion tip with your other hand while holding the lower part of the T part in the bag.
6. Insert the spike by ¾ of its length into outlet port. Do not touch the spike of the transfusion set and open outlet port to prevent contamination.
7. After inserting the spike, squeeze the bag to the level that transfusion set filter is filled with blood.
8. Hang the bag onto the stand from the bottom.

## GENERAL WARNINGS

### During storage and transportation:

1. During shipment, up to two blood bags shippers shall be carried at the same time manually or on a trolley. Parcels should not be dropped and/or crushed and should be loaded in a way to prevent fall over.
2. In the warehouse, up to eight parcels can be stored on top of each other and in a way to prevent fall over.

### Before opening the outer packaging:

1. Keep blood bags with overwrap at clean, dry and cool places.
2. Do not use if there is visible sign of deterioration on the overwraps.
3. Blood bags are packed with sterilized.

### After the outer packaging is opened:

1. Store the bags in the overwrap at cool and dry places.
2. Blood bags should be used within 10 days of opening their outer wraps.

### Before use:

1. Blood bags should be used in accordance with the instructions for use.
2. Intended for human blood and components.
3. Sterile and non-pyrogenic — Sterilized by Steam.
4. Do not use non-clear solution bags.
5. Do not use if there is any visible sign of deterioration or if any component is lose or missed.
6. Protect the bag and tubing from sharp objects.
7. Do not use if needle protector is opened. The lines on the needle should be directly opposite each other (Figure 11).
8. Do not close after opening the needle cover.
9. Do not use if there are no liquid route covers or if they are loose.
10. Do not aerate.
11. Store the blood between 1-6°C.
12. Do not add medicine to the blood bag.
13. Compare the blood in the bag and the recipient's blood with the Cross test before the transfusion.
14. Use a transfusion kit in accordance with EN ISO 1135-4.
15. Mix before use and, if necessary, bring to the appropriate temperature under controlled conditions.
16. When frozen, plastic is more fragile.


Figure 11

Figure 10

### After use

Single use only After usage it should be thrown as Medical Waste **CAUTION:** A risk of infection when used more than once

**CAUTION:** This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

KANSUK LABORATUARI SANAYİ VE TİCARET A.Ş.  
Yassıören Mahallesi Fırat Sk. 14/1 Arnavutköy-Istanbul/Turkey  
Phone: +90 212 592 15 76, Fax: +90 212 580 37 72  
www.kansuk.com E-mail: kansuk@kansuk.com

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