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
www.kansuk.com, E-mail: kansuk@kansuk.com

APPLICABLE CODES

20001046555	LABORATORY LEUKOREDUCTION FILTER FOR PLATELET, BUFFY COAT POOLING/ OCT
20004046555	LABORATORY LEUKOREDUCTION FILTER FOR PLATELET, BUFFY COAT POOLING/ OCT TYPE I
20004146519	LABORATORY LEUKOREDUCTION FILTER FOR PLATELET, BUFFY COAT POOLING/ TRAIN
20004146520	LABORATORY LEUKOREDUCTION FILTER FOR PLATELET, BUFFY COAT

DESCRIPTION

Octopus System

REF 20004046555 - 20001046555

This set contains one 600 mL pooling bag with 6 tubing leads for connection of buffy coats (BC) and one additional tubing lead for Platelet Additive Solution (PAS) dilution, an inline Sepacell™ PLX-5 platelet leukodepletion filter connected to 1,3 liter platelet storage bag with a 100 mL sampling pouch (Figure 1).


Figure-1

Train Systems

REF 20004146519 - 20004146520

This set contains one 600 mL pooling bag with one additional tubing lead for PAS dilution, an inline Sepacell™ PLX-5 platelet leukodepletion filter connected to 1,3 liter platelet storage bag with a 100 mL sampling pouch (Figure 2).

This set contains a tubing lead with inline Sepacell™ PLX-5 platelet leukodepletion filter connected to a 1,3 liter platelet storage bag with a 100 mL sampling pouch (Figure 3).


Figure-2


Figure-3

- As optional, in some systems there is pH sensor to be connected to the pH SAFE reader device to read pH of the platelets in the storage bag.

INTENDED USE

The Platelet Pooling Set with Leukocyte Filter is intended for the preparation of leukocyte depleted platelets obtained by pooling BCs derived from whole blood by adding plasma or additive solution (PAS). Validations were conducted using BC platelets prepared from pools of up to 4-6 BC.

PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by certified Blood Bank personnel in the Blood Bank

INDICATIONS FOR USE

This product is intended to be used for producing leukocyte-reduced platelets for all population to aid any treatment. Platelets can be used for restoring blood components and reducing and preventing bleeding according to the treatment decision of physician.

CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

PRECAUTIONS

- Before preparing BCs from whole blood, a strict standardized and validated procedure for the preparation of venipuncture site should be in operation including skin disinfection and blood diversion of the first aliquots in order to reduce the risk of bacterial contamination.
- Follow validated blood component preparation procedures per local requirements.
- Use transfusion sets in compliance with EN ISO 1135-4.
- The tubes are suitable to be connected with a Sterile Connecting Device.
- It is intended for human blood and components.
- During shipment, up to two containers may be carried at the same time manually or on a trolley, they must not be dropped or any heavy items loaded on top of them, and, they should be loaded in such a way to prevent them falling over.
- Single use only, after usage it should be thrown as a Medical Waste. **Caution:** Risk of infection when used more than once.
- For the explanation of symbols, see Symbols Description Table.

Caution: This medical device contains di(2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer

to current literature to make an informed decision.

PROCESS STEPS

A. GENERAL PREPARATION BEFORE POOLING

1. For optimizing the platelet yield in the platelet concentrate, pre-selection of BC is recommended.
2. Read manufacturer instruction for use for blood bags, centrifuge, blood component separator and procedures meeting local validation requirements to produce BC.
3. The recommended separation time for optimal filtration of the platelet concentrate is approximately 6 minutes. The blood component separator should be set up properly before use.
4. BC platelets should be prepared within 24 hours following whole blood collection.
5. Open the primary package and check the The Platelet Pooling Set with Leukocyte Filter for any visible defects. Do not use defective sets. **Caution:** The bag system may give off a different smell because of the Gamma sterilization.

B. POOLING OF BUFFY COATS

Octopus method of pooling (Recommended, Figure 1)

1. Sterile connect BC (4-6) to each dedicated tubing lead of the Octopus.
2. Close all the clamps.
3. Sterile connect a PAS/Plasma bag to the last tubing lead of the Octopus.
4. Hang the set.
5. Open the clamps and let all the BC flow into the pooling bag.
6. Squeeze the pooling bag gently, remove the air collected in the bag to the empty BC bag and close the clamp.
7. To maximize the BC content recovery, rinse each BC bag sequentially using a portion of the PAS fluid or Plasma. Let the remaining volume of PAS/Plasma flow into the pooling bag. Seal and detach the octopus from the pooling bag. **Caution:** The pooling bag should be filled with a maximum of 600 ml of BC and PAS/Plasma. Make sure that, the outlet ports are free from red blood cells (RCs) before centrifugation. RCs may pass into the platelet storage container or block the filter.

Train method of pooling (Recommended, Figure 2)

1. Sterile connect BC one to each other using the top and bottom tubings of the bags like a train.
2. Sterile connect the first BC of the train to a PAS bag or plasma.
3. Sterile connect the last (bottom) BC bag to the tube of the pooling bag of the pooling set leaving out 15-20 cm of tubing and close the clamp.
4. Place a re-usable clamp on the tube between the PAS/Plasma bag and the first BC bag.
5. Hang the system from the PAS/Plasma bag.
6. Let all BC flow by gravity into the pooling bag.
7. Squeeze the pooling bag gently, remove the air collected in the bag to the empty BC bag and close the clamp.
8. To maximize the collection of BC content, rinse each BC bag sequentially using a portion of the PAS fluid or plasma. Let the remaining volume of PAS/Plasma flow into the pooling bag. Seal and detach the rest of the train from the last bag containing the pool. **Caution:** Make sure that outlet ports are free from RCs before centrifugation. RCs may pass into the platelet storage container or block the filter.

Train method of pooling (Recommended, Figure 3)

1. Sterile connect BC one to each other using the top and bottom tubings of the bags like a train.
2. Sterile connect the first BC of the train to a PAS bag or plasma.
3. In order to use the last (bottom) BC bag as a pooling bag, sterile connect it to the tube of the pooling set leaving 15-20 cm of tubing and close the clamp.
4. Place a re-usable clamp on the tube between the PAS/Plasma bag and the first BC bag.
5. Hang the system from the PAS/Plasma bag.
6. Let all BC flow by gravity into the last bag of the train.
7. Squeeze the last pooling bag gently, remove the air collected in the bag to the empty BC bag and close the clamp.
8. To maximize the collection of BC content, rinse each BC bag sequentially using a portion of the PAS fluid or Plasma. Let the remaining volume of PAS/Plasma flow into the last bag of the train. Seal and detach the rest of the train from the last bag containing the pool. **Caution:** Make sure that outlet ports are free from RCs before centrifugation. RCs may pass into the platelet storage container or block the filter.

C. CENTRIFUGATION

1. Place the filter protective apparatus (Kansheet) onto the storage bag, fold them in the manner that the filter is covered and pack the assembly of pooled BCs into centrifuge buckets according to blood center validated procedure.
2. Choose the appropriate platelet separation program. Follow user instructions of the centrifuge manufacturer. **Caution:** Pack the sets into the centrifuge buckets carefully in order to avoid damages to the bags or poor separation results.

D. COMPONENT SEPARATION

1. After centrifugation, hang the set onto the pins of the blood component separator device.
2. Ensure the leukodepletion filter, tubings and platelet storage bag are placed on the device properly. Close the slide clamp of the sampling bag.
3. Break the breakaway part of the pooling bag by turning it left and right. Start the device using the appropriate platelet separation program. Follow user instructions of the manufacturer.
4. The press plate will press the collection bag and platelets will flow into the storage bag.
5. When RCs reach the tube downstream of the filter, stop the separation.
6. Seal and detach the platelet storage bag from the rest of the assembly.
7. To remove the air, gently fold the storage bag and squeeze the air into the sampling pouch.

E. STORAGE

1. Platelets can be stored at 22 °C under agitation in an incubator up to 5 days in PLT bag made of TOTM.
2. In containers suitable for storage up to 7 days, platelets can be stored up to 7 days provided that regular invasive controls of bacterial contamination are performed.
3. In containers equipped with pHsensor, by performing regular non-invasive controls of pH, the status of platelets can be determined. **Caution:** Please read BCSI pH1000 User manual instructions.

F. TRANSFUSION

1. Inspect the container and/or blood component for any defects.
2. Before transfusing the blood, make sure that serological test results are acceptable, cross match test results match with the recipient. If appropriate, apply cross match again at bedside.
3. Platelets should be inspected before transfusion to ensure that no platelet aggregates are visible.
4. Hold the top part of the T part with one hand and bag from the bottom part of T part with other hand, and twist of the upper part by 270°.
5. Keep holding the bag from T part, remove the cap of the transfusion spike with other hand.
6. Insert the spike by ¾ of its length into outlet port. Do not touch the spike of the transfusion set and open outlet port to prevent contamination.
7. After inserting the spike, squeeze the bag to the level that transfusion set filter is filled with blood.
8. Hook the bag from its holder. The system is ready for transfusion.

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