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Şekil 3

APPLICABLE CODES:

20000006517	BLOOD TRANSFER BAG FOR RED BLOOD CELL WASHING WITH NACL SOLUTION
-------------	--

DESCRIPTION

A blood bag system made of PVC and connected in a closed system with transfer bags including NaCl solution for washing.


Figure 1

INTENDED USE

It is used for washing erythrocyte suspension with serum.

PATIENT POPULATION

The device is used in paediatric and adult populations for blood and blood component transfusion on the order of a physician. Blood can be donated by adults who fulfill national/international health criteria.

INTENDED USER AND INTENDED USE ENVIRONMENT

This product is intended to be used by certified Blood Bank personnel in the Blood Bank.

INDICATIONS FOR USE

This product is intended to be used for washing of red blood cells for all population to aid any treatment. Red blood cells can be used for restoring blood components and improving oxygen delivery to tissues according to the treatment decision of physician.

CONTRAINDICATIONS

This product is not intended to be reused or re-sterilized. The risks that may arise from reuse or re-sterilization (e.g. infection, cross contamination, loss of performance) could result in serious injury or death.

GENERAL PREPARATIONS BEFORE USAGE


BLOOD TRANSFER TO THE BAG BEFORE WASHING

1. Before washing, put the cell washer kit in the + 4 / 6 °C blood storage container and bring the temperature of the serum to the same temperature as the blood component to be washed.
2. Remove the blood bag from its outer package.
3. Before starting blood transfer operation, inspect the bag system, tubings and plastic spike in terms of visible errors and color changes. Apply pressure to detect any leakage. Check that the NaCl solution is filled in the correct volume and does not contain color or particles.
4. Identify the bag to be used for transferring erythrocyte suspension to be washed.

CAUTION: The tube of the bag in which the blood to be washed is located must be long.

5. The plastic spike on the bag to which the blood is to be transferred is closed by tube closure and released from the system.
6. Washing erythrocyte suspension / tube of whole blood bag and tube of transfer bag to quadruple cell washing is placed in the device.
7. According to the instructions for use of sterile hose assembling device, two bag tubes are assembled in sterile environment.
8. The bag which is connected sterile to the system and has blood to be washed is hung and the blood is transferred to the quadruple cell washing system via the gravity.
9. The empty bag is closed and released by the tube closing device from the system after the transfer process.

Figure 2

INSTRUCTIONS FOR USE

BLOOD WASHING

1. Break the breakaway part of the first bag (Figure 1, No: 1) which contains red blood cell suspension and the saline bag (Figure 1, No: 2) which contains NaCl 0.9% as washing solution (Figure 2) and add 250 ml of the solution onto the blood and gently mix.
2. After transfer, the tube of the bag containing 0.9% NaCl solution is closed with clamp.
3. The erythrocytes transferred to the 0.9% NaCl solution are centrifuged at 1200 G, Duration: 10 minutes, Acceleration: 9, Slow: 4, Temperature: + 1-6 °C.
CAUTION: Centrifuge settings can vary according to the centrifuge brand.
4. In the first centrifuge bag, the light red solution formed on the top is transferred to the empty transfer bag (Figure 1, No: 3).
5. After this process, the bag containing the light red solution is closed with the tube closing device from the system and separated.
6. Open the clamp of bag again the 0.9% NaCl solution on the erythrocytes and add the remaining 250 ml of 0.9% NaCl solution and mix gently.
7. After centrifugation, the light red solution formed on the top of the first bag is transferred to the empty transfer bag (Figure 1, No: 4).
8. After this process, the bag containing the light red solution is closed with the tube closing device from the system and separated.
9. Finally, one unit of washed RBC solution is ready to be used.
10. Erythrocytes that have been washed with 0.9% NaCl solution should be used within 24 hours of washing.

TRANSFUSION

1. Inspect the container and blood component for any defects.
2. Before transfusing the blood, make sure that serological test results are acceptable, cross match test results match with the recipient. If appropriate, apply cross match again at bed side.

3. Use transfusion set compatible with EN ISO 1135-4.
4. Hold the top part of the T part with one hand and bag from the bottom part of T part with other hand and twist of the upper part by 270° (Figure 3).
5. Keep holding the bag from T part, remove the cap of the transfusion spike with other hand.
6. Insert the spike by ¾ of its length into outlet port. Do not touch the spike of the transfusion set and open outlet port to prevent contamination.
7. After inserting the spike, squeeze the bag to the level that transfusion set filter is filled with blood.
8. Hook up the bag from its holder.

GENERAL CAUTIONS

During storage and transportation:

1. During shipment, up to two blood bags shippers shall be carried at the same time manually or on a trolley. Parcels should not be dropped and / or crushed and should be loaded in a way to prevent fall over.
2. The products should be stored in their original unopened / undamaged packaging at 0- +25 ° C, protected from direct sunlight.
3. In the warehouse, up to eight parcels can be stored on top of each other and in a way to prevent fall over.

Before opening the overwrap:

3. Keep blood bags with overwrap at clean, dry and cool places.
4. Do not use if there is visible sign of deterioration on the overwraps.
5. Blood bags are packed with sterilized.

After opening overwrap:

6. Store the bags in the overwrap at cool and dry places.
7. Blood bags should be used within 10 days of opening their outer wraps.

Before usage:

8. Please follow the user instructions.
9. Intended for human blood and blood components.
10. Sterile and non-pyrogenic – Sterilized by Steam.
11. Do not use non-clear solution bags.
12. Do not use if there is any visible sign of deterioration or if any component is lose or missed.
13. Protect the bag and tubing from sharp objects.
14. Do not use if the fluid path lids are not available or not tight.
15. Do not give any air.
16. Do not add medicine into the blood bag.
17. Compare the pre-transfused blood and the recipient's blood with the cross test.
18. Use transfusion set compatible with EN ISO 1135-4.
19. Stir before use and, if necessary, bring to appropriate temperature under controlled conditions.
20. When frozen, plastic is more fragile.

After usage:

Single use only. After usage it should be thrown as Medical Waste. **CAUTION:** A risk of infection when used more than once.

CAUTION: This medical device contains di (2-ethylhexyl) phthalate (DEHP). Children including developing foetus, premature neonates and neonates and therefore potentially pregnant and nursing women are at higher risk to DEHP exposure, although humans have been observed to be less sensitive than rodents. Furthermore, exchange transfusion of blood in neonates, extracorporeal membrane oxygenation treatment of neonates and adults, total parenteral nutrition in neonates, haemodialysis, enteral nutrition in neonates and adults, heart transplantation or coronary artery bypass graft surgery, massive blood transfusion of red blood cells and plasma or peritoneal dialysis may lead to high exposure to DEHP. Medical procedures should not be avoided however, as the benefits outweigh any possible health risks associated with DEHP exposure. Please refer to current literature to make an informed decision.

Figure 3

KANSUK LABORATUVARI SANAYİ VE TİCARET A.Ş.
Yassören Mahallesi Fırat Sk. No:14/1 Arnavutköy-İstanbul/Turkey
Tel: +90 212 592 15 76, Faks: +90 212 580 37 72,
www.kansuk.com, e-mail: kansuk@kansuk.com

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