


ACD-A ANTICOAGULANT SOLUTION 500ML



| | |
|---------------------------------|--|
| Products | 20000006525 ACD-A ANTICOAGULANT SOLUTION 500ML |
| Sterilization Method | Steam Sterilization |
| Shelf Life | 24 Months |
| Packaging Form | Primary Packaging (Transparent): 1 unit 20 units in 1 box |
| Product Description | 500 ml ACD-A anticoagulant solution contained in transparent PVC bags within a transparent package. |
| Intended Use /Indication | According to the "Guideline for the Preparation, Use and Quality Assurance of Blood Components", the ACD-A solution is intended to be used as an anticoagulant solution during the collection of human blood in apheresis devices. |
| Performance | The anticoagulant solution ACD-A prevents coagulation by binding to active calcium ions through the citrate it contains. |

KT.038.2026.00- 26.03.2026



FOR DETAILED INFORMATION, PLEASE CONTACT

Kansuk Laboratuvarı Sanayi ve Ticaret A.Ş. Yassiören Mahallesi Fırat Sk. No:14/1 Arnavutköy-İstanbul/Turkey

Phone: +90 212 5921576 Fax: +90 212 5803772 Web: www.kansuk.com E-mail: kansuk@kansuk.com

